

## Product datasheet for **MR227637A1V**

### Mouse Psma6 (NM\_011968) AAV Particle

#### Product data:

|                           |                                                                             |
|---------------------------|-----------------------------------------------------------------------------|
| Product Type:             | AAV Particles                                                               |
| Product Name:             | Mouse Psma6 (NM_011968) AAV Particle                                        |
| Tag:                      | Myc-DDK                                                                     |
| Symbol:                   | Psma6                                                                       |
| Synonyms:                 | IOTA                                                                        |
| Mammalian Cell Selection: | None                                                                        |
| Vector:                   | pAAV-AC-Myc-DDK (PS100089)                                                  |
| ORF Nucleotide Sequence:  | >MR227637 representing NM_011968<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGTCCCGGGTTCCAGCGCCGGTTTTGACCGTCACATTACTATCTTCTCTCCCAGGGCCGACTCTACC  
AAGTAGAATATGCTTTTAAGGCTATTAACCAGGTGGACTTACATCTGTTGCTGTCAGAGGAAAGGACTG  
TGCAGTTATTGTCACACAGAAGAAAGTACCCGACAACTACTGGATTCCAGCACAGTGACTCACTTGTTT  
AAGATAACGGAAAGCATTGGCTGTGTGATGACAGGAATGACAGCTGACAGCAGATCCCAGGTACAGAGGG  
CACGCTATGAGGCAGCTAATTGGAAATACAAGTATGGCTATGAAATCCTGTGGACATGCTGTGCAAAAG  
AATTGCTGATATTTCTCAAGTCTACACACAGAACGCTGAAATGAGGCCACTCGGTTGTTGTATGATTTTA  
ATCGGTATCGATGAAGAGCAAGGCCCTCAAGTGTACAAGTGTGATCCTGCAGGTTACTACTGTGGCTTTA  
AAGCCACGGCGGAGGAGTTAAACAAACAGAGTCAACCAGCTTCCTTGAAAAGAAAGTGAAGAAGAAATT  
TGATTGGACATTTGAACAGACAGTGAAACTGCAATTACATGCCTGTCTACTGTTCTGTCGATTGACTTC  
AAACCTTCAGAAATAGAAGTTGGAGTAGTTACAGTTGAAAATCCTAAATTCAGGATTCTTACTGAAGCAG  
AGATTGACGCTCACCTTGTTGCTCTAGCAGAGAGAGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Protein Sequence:</b> | <p>&gt;MR227637 representing NM_011968<br/> <span style="color: red;">Red</span>=Cloning site <span style="color: green;">Green</span>=Tags(s)</p> <p>MSRGSSAGFDRHITIFSPEGRLYQVEYAFKAINQGGLTSVAVRGKDCAVIVTQKKVPDKLLDSSTVTHLF<br/> KITESIGCVMTGMTADSRSQVQRARYEAAANWKYKYGYEIPVDMLCKRIADISQVYTQNAEMRPLGCCMIL<br/> IGIDEEQGPQVYKCDPAGYYCGFKATAAGVKQTESTSFLEKKVKKKFDWTFEQTVETAITCLSTVLSIDF<br/> KPSEIEVGVVTVENPKFRILTEAEIDAHLVALAERD</p> <p><span style="color: red;">TR</span><span style="color: green;">TRPLEQKLISEEDLAANDILDYKDDDDKV</span></p> |
| <b>Species:</b>          | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Serotype:</b>         | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>ACCN:</b>             | NM_011968                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>ORF Size:</b>         | 738 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Buffer:</b>           | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Stability:</b>        | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.                                                                    |
| <b>RefSeq:</b>           | <u>NM_011968.3, NP_036098.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq Size:</b>      | 997 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq ORF:</b>       | 741 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Locus ID:</b>         | 26443                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>UniProt ID:</b>       | <u>Q9QUM9</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>     | 12 C1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>MW:</b>               | 27.8 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |