

Product datasheet for **MR227677L1V**

Bmp2 (NM_007553) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bmp2 (NM_007553) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bmp2
Synonyms:	AI467020; Bmp; Bmp2a
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_007553
ORF Size:	1182 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227677).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007553.2
RefSeq Size:	2540 bp
RefSeq ORF:	1185 bp
Locus ID:	12156
UniProt ID:	P21274
Cytogenetics:	2 65.21 cM



[View online »](#)

Gene Summary:

This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer, which plays a role in bone and cartilage development. While a homozygous knockout mouse for this gene exhibits embryonic lethality, conditional knockout mice have defects in bone, cartilage and heart development. [provided by RefSeq, Jul 2016]