

## Product datasheet for **MR227679**

### **Cbl (NM\_007619) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Cbl (NM_007619) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cbl
Synonyms:	4732447J05Rik; c-Cbl; Cbl-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR227679 representing NM\_007619  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGCCGGCAACGTGAAGAAGAGCTCGGGCGCCGGCGGGCGGGCTCTGGGGCTCGGGAGCGGGCGGCC  
 TGATCGGGCTCATGAAGGACGCCTTCCAGCCGCACCACCACCACCACCCTCAGCCCGCACCTCCCTG  
 CACGGTGGACAAGAAGATGGTGGAGAAGTCTGGAAGCTCATGGACAAGGTGGTGCGGTTGTGTCAAAAC  
 CCAAAGCTGGCGCTCAAGAACAGCCCGCTTATATCTTAGACCTGCTGCCTGACACCTACCAGCACCTCC  
 GCACTGTCTTGTCAAGATATGAGGGGAAGATGGAGACGCTTGGAGAAAATGAGTATTTACAGGGTGTTCAT  
 GGAAAATTTGATGAAGAAAATAAGCAGACTATCAGCCTCTTCAAGGAGGGAAAAGAAAGGATGTATGAG  
 GAGAATCCCAGCCTAGGCGAAACCTGACCAAATTATCCCTGATCTTCAGCCACATGCTGGCAGAAGCTGA  
 AAGGCATCTTCCGAGCGGACTTCCAAGGAGACACTTCCGGATTACTAAAGCTGATGCTGCCGAATT  
 TTGGAGAAAAGCTTTTGGTGAAGAACGATAGTCCCGTGAAGAGCTTTCGACAGGCCCTGCATGAAGTG  
 CATCCCATCAGTTCTGGGCTGGAGGCCATGGCTCTGAAGTCCACTATTGATCTGACTGCAATGATTATA  
 TTTCTGTCTTTGAATTTGATATTTTACACGGCTGTTTACGCCCTGGTCCTCTTTGCTCAGAAAATTGGAA  
 CAGCCTTGCTGTAACCTACCCTGGTTACATGGCTTTCCTGACATACGATGAAGTGAAGCGCGCCTGCAG  
 AAGTTCATCCACAACCTGGCAGTTACATCTTTCGGCTGAGCTGTACTCGTTTGGGTGAGTGGGCTATTG  
 GGTATGTTACTGCCGATGGGAACATTCTGCAGACAATCCACACAATAAACCGCTCTTCCAAGCACTGAT  
 TGATGGCTTCAGGAAGGCTTCTATTTGTTTCTGATGGACGAAATCAAATCCTGACCTGACAGGTTTA  
 TGTGAACCAACTCCTCAAGATCATATCAAAGTAAACCCAGGAACAATATGAATTACTGTGAAATGGGCT  
 CCACATTTCAACTGTGAAGATATGTGCTGAGAATGATAAGGATGTGAAGATTGAGCCCTGTGGCCACT  
 CATGTGCACATCCTGCCTCACGTCGTGGCAGGAATCAGAAGGTCAGGGCTGTCCTTTTTGCCGATGTGAA  
 ATCAAAGTACTGAGCCATCGTGGTGGATCCGTTTGACCCAGAGGCAGTGGCAGCCTATTAAGGCAAG  
 GAGCAGAAGGTGCTCCTTCCCCAAATTACGACGATGATGATGATGAACGAGCTGATGATTCTCTTTCAT  
 GATGAAGGAGTTGGCAGGTGCCAAGGTGGAAGGCCTTCTCTCCATTCTCCATGGCCCCACAAGCTTCC  
 CTTCTCCAGTGCCACCAAGACTTGACCTTCTACAGCAGCGAGCACCTGTTCTGCCAGCACTTCAGTTC  
 TGGGGACTGCTTCCAAGGCTGCTTCTGGCTCCCTTATAAAGACAACCATTGCCAATACCTCCCACACT  
 TCGAGATCTTCCACCACCACCCCTCCAGACCGGCTTACTCTGTTGGAGCAGAAAACAAGGCCTCAGAGA  
 CGCCCTCTGCCTGTACACCAGGCGATTGTCCATCTAGAGACAACTGCCCCCTGTCCCTCTAGCCGCC  
 CAGGGGACTCGTGGTTGTCTCGGCAATCCCTAAAGTACCAAGTACTCCTCAAACCCCTGGTGATCCTTG  
 GAATGGGAGAGAATTGACCAATCGGCACTCGCTTCCATTTCTATTGCCCTCACAATGGAACCCAGAGCA  
 GATGTCCCTAGGCTTGGAAAGCACATTTAGTCTGGATACCTCTATGACTATGAATAGCAGCCAGTAGCAG  
 GTCCAGAGAGTGAGCACCCAAAGATCAAGCCTTCTCTGCTGCCAACGCCATTTACTCTCTGGCTGCCAG  
 GCCTCTTCTATGCCAAAAGTCCACCTGGGGAGCAAGGGGAAAGTGAAGAGGACACAGAATATATGACT  
 CCCACATCTAGGCCTGTAGGGGTTGAGAAGCCAGAGCCCAACCGCCGTTAGAGGCAACCCAGAGTTCAC  
 GAGCATGTGACTGTGACCAGCAGATCGACAGCTGTACCTATGAAGCGATGTATAACATCCAGTCCAAGC  
 ACTGTCTGTAGCAGAAAACAGCGCCTCTGGGGAAGGGAATCTGGCCACAGCTCACACGAGTACTGGCCCT  
 GAGGAATCCGAAAACGAGGATGATGGCTATGATGTGCCTAAGCCACCCGTCAGCTGTACTGGCCCGCC  
 GGACCCTGTCTGACATCTCCAATGCCAGCTCCTCCTTTGGCTGGTTGTCTTTGGATGGTGACCCTACAAA  
 CTTCAATGAGGGTTCCCAAGTCTCTGAGCGGCCCCCAACCAATTCCCTCGGAGAATCAACTCAGAACGA  
 AAAGCCAGTAGCTATCAACAAGGCGGAGGTGCCACTGCTAACCCGTGGCCACAGCACCCCTACCAGCAGC  
 TCTCAAGTGAATTGAACGCCTCATGAGTCAGGGCTATTCTACCAGGACATTAGAAAAGCTTTGGTTCAT  
 TGCCCAACAACATTGAGATGGCTAAAAACATCCTCCGGGAATTTGTTTCTATTCTCTCTCTGCTCAC  
 GTAGCCACC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAAGTTTAA

**Protein Sequence:** >MR227679 representing NM\_007619  
 Red=Cloning site Green=Tags(s)

```
MAGNVKKSSGAGGGGSGGSGAGGLIGLMKDAFQPHHHHHHLSHPHPCTVDKMKVEKCWKLMDKVVRLCQN
PKLALKNSPPYILDLLPDTYQHLRTVLSRYEGKMETLGENEYFRVFMENLMKKTKQTIISLFKEGKERMIE
ENSQPRRNLTKLSLIFSHMLAELKGIFFSGLFQGDTRITKADAAEFWRKAFGEKTIIVPWKSRQALHEV
HPISSGLEAMALKSTIDLTCNDYISVFEFDIFTRLFQPWSSLLRNWNSLAVTHPGYMAFLTYDEVKARLQ
KFIHKPGSYIFRLSCTRLGQWAIQVYVADGNILQTIHPNKPLFQALIDGFREGFYLPDGRNQNPDLTGL
CEPTPQDHIKVTQEYELYCEMGSTFQLCKICAENDKDKVIEPCGHLMCTSCLTSWQESEQGCPFCRCE
IKGTEPIVDPDFPRGSGLLRQGAEGAPSPNYDDDDERADDSLFMMKELAGAKVERPSSPFSMAPQAS
LPPVPPRLDLLQQRAPVPASTSVLGTASKAASGSLHKDKPLPIPPTLRDLPPPPPPDRPYSVGAETRPQR
RPLPCTPGDCPSRDKLPVPSSRPGDSWLSRPIKVPVATPNPGDPWNGREL TNRHSLPFSLPSQMEPRA
DVPRLGSTFSLDTSMTMNSSPVAGPESEHPKIKPSSANAIYSLAARPLPMPKLPPGEQGESEEDTEYMT
PTSRPVGVQKPEPKRPLEATQSSRACDCDQQIDISCTYEAMYNIQSQAL SVAENSASGEGNLATAHTSTGP
EESENEDDGYDVPKPPVAVLARRTLSDISNASSFGWLSLDGDPNTFNEGSQVPERPPKPFPRRINSER
KASSYQQGGGATANPVATAPSPQLSSEIERLMSQGYSDIQKALVIAHNNIEMAKNILREFVSISSPAH
VAT
```


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9032\\_e12.zip](https://cdn.origene.com/chromatograms/mm9032_e12.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:




\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_007619

**ORF Size:** 2739 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_007619.2</a> , <a href="#">NP_031645.2</a>
<b>RefSeq Size:</b>	5083 bp
<b>RefSeq ORF:</b>	2742 bp
<b>Locus ID:</b>	12402
<b>UniProt ID:</b>	<a href="#">P22682</a>
<b>Cytogenetics:</b>	9 24.72 cM
<b>MW:</b>	101 kDa
<b>Gene Summary:</b>	Adapter protein that functions as a negative regulator of many signaling pathways that are triggered by activation of cell surface receptors. Acts as an E3 ubiquitin-protein ligase, which accepts ubiquitin from specific E2 ubiquitin-conjugating enzymes, and then transfers it to substrates promoting their degradation by the proteasome. Recognizes activated receptor tyrosine kinases, including KIT, FLT1, FGFR1, FGFR2, PDGFRA, PDGFRB, EGFR, CSF1R, EPHA8 and KDR and terminates signaling. Recognizes membrane-bound HCK, SRC and other kinases of the SRC family and mediates their ubiquitination and degradation. Participates in signal transduction in hematopoietic cells. Plays an important role in the regulation of osteoblast differentiation and apoptosis. Essential for osteoclastic bone resorption. The 'Tyr-737' phosphorylated form induces the activation and recruitment of phosphatidylinositol 3-kinase to the cell membrane in a signaling pathway that is critical for osteoclast function. May be functionally coupled with the E2 ubiquitin-protein ligase UB2D3 (PubMed:10393178, PubMed:12649282, PubMed:19265199, PubMed:20100865, PubMed:9653117). In association with CBLB, required for proper feedback inhibition of ciliary platelet-derived growth factor receptor-alpha (PDGFRA) signaling pathway via ubiquitination and internalization of PDGFRA (PubMed:29237719).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR227679