

Product datasheet for **MR227712**

Nrep (NM_001267717) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Nrep (NM_001267717) Mouse Tagged ORF Clone

Tag: Myc-DDK

Symbol: Nrep

Synonyms: AI325076; D0H4S114; Harp; P311; PTZ17; SEZ17

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR227712 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTTTACTACCCAGAACTCTTGGTCTGGTCAGCCAAGAACCCTTGCATACAAGGAAATGGAGGGAG
GTCTTATTAAGGGAAGACTTCCCGTGCCTAAGGAAGTGAACCGAAAGAAGATGGAGGAGACTGGCGCTGC
CTCGTGACTCCGCCAGGCAGCCGTAATTCACCTCTCCAGCTACCAGTTACCTCCACCCCTTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227712 protein sequence
Red=Cloning site Green=Tags(s)

MVYYPELLVWVSQEPFAYKEMEGGLIKRLLPVPKEVNRKKMEETGAASLTPPGSREFTSPATSYLHPF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001267717

ORF Size: 204 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001267717.1](#), [NP_001254646.1](#)

RefSeq Size: 2097 bp

RefSeq ORF: 207 bp
 Locus ID: 27528
 UniProt ID: [Q07475](#)
 Cytogenetics: 18 B1
 MW: 7.7 kDa

Gene Summary: May have roles in cellular differentiation. Ectopic expression induces differentiation of fibroblast into myofibroblast and myofibroblast ameboid migration. Increases retinoic-acid regulation of lipid-droplet biogenesis. May also have neural functions. Promotes axonal regeneration and augments motility of gliomas. Down-regulates the expression of TGFB1 and TGFB2 but not of TGFB3. May play a role in the regulation of alveolar generation. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR227712