

## Product datasheet for **MR227718A1V**

### Mouse Gngt2 (NM\_001284393) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Gngt2 (NM\_001284393) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Gngt2  
**Synonyms:** AV096488  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR227718 representing NM\_001284393  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

**ATG**CCCCAGGACCTCAGTGAGAAGGAGCTGTTGAGGATGGAGGTGGAGCAGCTGAAGAAGGAAGTGAAGA  
ACCCACGTGATCTGATTTCCAAGACAGGAAAGGAAATCAAGGATTATGTAGAGGCCCAAGCAGGGACAGA  
CCCTCTTCTCAAAGGCATCCAGAAGACAAGAATCCCTTCAAGGAGAAAGGCACCTGTGTGCTAAGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR227718 representing NM\_001284393  
**Red**=Cloning site **Green**=Tags(s)

MAQDLSEKELLRMEVEQLKKEVKNPRDLISKTKKEIKDYVEAQAGTDPLLKGIPEKPNPFKEKGTCVLS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_001284393  
**ORF Size:** 207 bp  
**Buffer:** PBS with 0.001% Pluronic F68



[View online »](#)

**Stability:** AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

**RefSeq:** [NM\\_001284393.1](#), [NP\\_001271322.1](#)

**RefSeq Size:** 590 bp

**RefSeq ORF:** 210 bp

**Locus ID:** 14710

**UniProt ID:** [Q61017](#)

**Cytogenetics:** 11 59.01 cM

**MW:** 7.8 kDa