

Product datasheet for **MR227721A1V**

Mouse Gng7 (NM_010319) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Gng7 (NM_010319) AAV Particle
Tag: Myc-DDK
Symbol: Gng7
Synonyms: AI840417
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227721 representing NM_010319
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGTCAGGTACTAACAACGTCGCCAGGCCCGGAAGCTGGTGGAGCAGCTGCGCATTGAAGCTGGGA
TCGAACGCATCAAGGCTCCAAGGCCTCGTCAGACCTGATGGGCTACTGTGAGCAACATGCCCGCAACGA
TCCCTTGCTGGTTGGCGTGCCAGCCTCTGAGAATCCATTCAAAGACAAAAGCCTTGATAATTCTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227721 representing NM_010319
Red=Cloning site Green=Tags(s)
MMSGTNNVAQARKLVEQLRIEAGIERIKVSKASSDLMGYCEQHARNDPLLVGVPASENPFKDKKPCIIIL
TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_010319
ORF Size: 204 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_010319.3](#)

RefSeq Size: 3362 bp

RefSeq ORF: 210 bp

Locus ID: 14708

Cytogenetics: 10 39.72 cM

MW: 7.5 kDa