

Product datasheet for **MR227727**

Lsm2 (NM_001204274) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lsm2 (NM_001204274) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lsm2
Synonyms:	D17H6S56E-2; D17H6S56E2; Dmapl; Dmpkap; G7b; Sm-X5; SmX5; snRNP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR227727 representing NM_001204274 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC CGATCGCC ATGCTCTTCTACTCCTTTTTCAAGTCTCTTGTGGCAAGGATGTAGTCGTGGAAGTCAAGAATGACCTGA GCATCTGTGGAACCTCCACTCTGTGGACCAGTACCTCAATATCAAATAACCGACATCAGTGTCACAGA CCCTGAGAAATACCCTCACTTCCCATTGGTGACCCAGAAGTCTGACGCTCTACTGAGGACAAGAACCC CTGCCAATACT ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>MR227727 representing NM_001204274 Red=Cloning site Green=Tags(s) MLFYSFFKSLVGKDVVVELKNDLSICGTLHSVDQYLNKLTDISVTDPEKYPHFPIGDPELDVLLRTRTP LPNT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Restriction Sites:	Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001204274

ORF Size: 222 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001204274.1</u> , <u>NP_001191203.1</u>
RefSeq Size:	713 bp
RefSeq ORF:	225 bp
Locus ID:	27756
Cytogenetics:	17 18.52 cM
MW:	8.9 kDa
Gene Summary:	Plays role in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex that is involved in spliceosome assembly, and as component of the precatalytic spliceosome (spliceosome B complex). The heptameric LSM2-8 complex binds specifically to the 3'-terminal U-tract of U6 snRNA.[UniProtKB/Swiss-Prot Function]