

Product datasheet for **MR227760A1V**

Mouse Anapc11 (NM_001038230) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Anapc11 (NM_001038230) AAV Particle
Tag: Myc-DDK
Symbol: Anapc11
Synonyms: 1110011I19Rik; R75218
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227760 representing NM_001038230
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGGTGAAAATTAATGTTGGAATGGTGTGGCCACTTGGCTCTGGGTAGCCAATGATGAGAACTGCC
GCATCTGCAGGATGGCGTTTAAATGGCTGCTGTCCAGACTGTAAGGTGCCCTGGTGATGACTGCCCTCGT
GTGGGGACAGTGCTCCCACTGCTTCCACATGCACTGCATCCTCAAGTGGCTGAATGCCGACGAGGTGCAG
CAGCACTGCCCATGTGTCGCCAGGAGTGAAGTTCAAAGAG

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227760 representing NM_001038230
Red=Cloning site Green=Tags(s)

MKVKIKCWNGVATWLWVANDENCGICRMAFNCCPDCKVPGDDCPLVWGQCShCFHMHCIKWLNAQQVQ
QHCPMCRQEWKFKE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001038230
ORF Size: 255 bp
Buffer: PBS with 0.001% Pluronic F68



[View online >](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001038230.2</u> , <u>NP_001033319.1</u>
RefSeq Size:	3384 bp
RefSeq ORF:	255 bp
Locus ID:	66156
UniProt ID:	<u>Q9CPX9</u>
Cytogenetics:	11 E2
MW:	9.8 kDa