

## Product datasheet for **MR227783**

### Banf1 (NM\_001286608) Mouse Tagged ORF Clone

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | Expression Plasmids   |
| Product Name:             | Banf1 (NM_001286608) Mouse Tagged ORF Clone                       |
| Tag:                      | Myc-DDK   |
| Symbol:                   | Banf1   |
| Synonyms:                 | Baf; Bcrp1; L2bp1/Baf   |
| Mammalian Cell Selection: | Neomycin  |
| Vector:                   | pCMV6-Entry (PS100001)  |
| E. coli Selection:        | Kanamycin (25 ug/mL)  |
| ORF Nucleotide Sequence:  | >MR227783 ORF sequence<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACAACCTCCCAAAGCACCGAGACTTCGTGGCAGAGCCCATGGGGAAAAGCCAGTGGGGAGCCTGG  
CCGGGATTGGTGACGTCTGAGCAAGAGGCTGGAGGAAAGGGCTTTGACAAGGCTTATGTGGTCTTGG  
CCAGTTTCTGGTGCTAAAGAAAGATGAAGACCTTCCGAGAATGGCTGAAGGATACATGTGGTCCCAAT  
GCCAAGCAGTCCCGGACTGCTTTGGTGCCTTCGAGAATGGTGTGATGCCTTCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

|                   |  |
|-------------------|--|
| Protein Sequence: | >MR227783 protein sequence<br>Red=Cloning site Green=Tags(s)                                     |
|                   | MTTSQKHRDFVAEPMGEKPVGSLAGIGDVL SKRLEERGFDKAYVVLGQFLVLLKDEDLFWREWLKDTCGAN<br>AKQSRDCFGCLREWCD AFL |

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

|                    |           |
|--------------------|-----------|
| Restriction Sites: | Sgfl-MluI |
|--------------------|-----------|



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001286608

**ORF Size:** 270 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001286608.1](#), [NP\\_001273537.1](#)

**RefSeq Size:** 1108 bp

**RefSeq ORF:** 270 bp

Locus ID: 23825


UniProt ID: [O54962](#)

Cytogenetics: 19 A

MW: 10.1 kDa

**Gene Summary:** Plays fundamental roles in nuclear assembly, chromatin organization, gene expression and gonad development. May potentially compress chromatin structure and be involved in membrane recruitment and chromatin decondensation during nuclear assembly. Contains 2 non-specific dsDNA-binding sites which may promote DNA cross-bridging.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR227783