

Product datasheet for **MR227834A1V**

Mouse Ncmap (NM_001243305) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ncmap (NM_001243305) AAV Particle
Tag: Myc-DDK
Symbol: Ncmap
Synonyms: A330049M08Rik; AA545175; Mp11
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227834 representing NM_001243305
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

**ATGACCACAGCCACCACTTTGGGGGACGCCGTCTTCTCGCTGAACATGACCAGGGGAGAGGACGCCCTGT
ATAAGAGCTCTGGAGCCATCGTGGCTGCCATTGTGGTGGTTGTTCATCATCATTGTACCTTGGTCTGAT
CCTGCTGAAGATGTACAACAGGAGAATGAGGACCAGGCCGGGAGCTGGAGCCCAAGAGCCCAAAGCCACCT
GTCCCTCTGCCCTGGACCCAGCAGCAATGGCAGCCAGCAGCCTGCTACTGTGACCTTCGACCTGCCA
ATGTCCATGTGGAGACTCGG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227834 representing NM_001243305
Red=Cloning site **Green**=Tags(s)

MTTATTLGDAVFLNMTRGEDALYKSSGAIVAIVVVVVIIVTLVLILLKMYNRRMRTRRELEPKSPKPP
VPPALDPSSNGSQPATVTFDPANVHVETR

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001243305
ORF Size: 300 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001243305.1</u> , <u>NP_001230234.1</u>
RefSeq Size:	1448 bp
RefSeq ORF:	303 bp
Locus ID:	230822
UniProt ID:	<u>Q99JS0</u>
Cytogenetics:	4 D3
MW:	10.8 kDa