

Product datasheet for **MR227859A1V**

Mouse H2a1f (NM_001242951) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse H2a1f (NM_001242951) AAV Particle
Tag: Myc-DDK
Symbol: H2a1f
Synonyms: Gm14479; OTTMUSG00000016785
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227859 representing NM_001242951
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCAAGAAAATGCAAAGGCGAAGAAGACAGAAAAGAACCCGCTCCCAGAGAGGTGAGCTTCCTTTCA
GCCTTGATAGTCGTTTCTACGAGAGGAATCCATTCCAGTCGCTGAGCTCTCCGCACTTTCATTCT
CAGGAGTGTGCTCGAGTACTTAACATCGAACATCCTCGAACTGGCTGGTGAGGTGGCCAGACCCTGGC
AGGAAGCGCATAGCTCCAGAGGATGTACGTCTGGTGGTACAGAACAACGAACAGCTCCGCCAACTCTTCA
AACCAGGTGGCACATCAGTGAATGAGGATGACAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227859 representing NM_001242951
Red=Cloning site Green=Tags(s)

MAKKMQRRRRQKRTRSQRGELPFSLVDRFLREEFHSSRLSSSALSFLTSTVLEYLTSNILELAGEVAQTTG
RKRIAPEDVRLVVQNNEQLRQLFKPGGTSVNEDDN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001242951
ORF Size: 318 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001242951.1, NP_001229880.1</u>
RefSeq Size:	318 bp
RefSeq ORF:	318 bp
Locus ID:	100042939
UniProt ID:	<u>Q5M8Q2</u>
Cytogenetics:	X A1.1
MW:	12.1 kDa