

Product datasheet for **MR227941A1V**

Mouse Gm2506 (NM_001199965) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Gm2506 (NM_001199965) AAV Particle
Tag: Myc-DDK
Symbol: Gm2506
Synonyms: 100039939
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227941 representing NM_001199965
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGATGGAGGGGCTCTCCCCGCCAGCAGCCTCCCGCTGTTACTGTTGCTTCTGAGCCCGGCTCCTGAAG
CAGCCTTGCTCTGCCCTCCAGCACTAGCTGCTGACTCAGCTCTATAGACAGCCACTCCAAGCAGGCT
GCTGAGGAGGATTGTCCACATGGAATCGAGGAGGCTGATGGGGACTGTACCTCCAGGCTGTCGTGCTT
CACCTGGCTCGGCGCAGTGTCTGTGTTCCATCCCAGAACCGCAGCCTGGCTCGGTGGTTAGAACGCCAAG
GGAAAAGGCTCCAAGGAAGTGTACCCAGTTTAAATCTGGTACTACAAAAGAAAATGTACTCACACCCCA
ACAGCAAAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227941 representing NM_001199965
Red=Cloning site Green=Tags(s)

MMEGLSPASSLPLLLLLLSPAPEAALPLPSSTSCCTQLYRQPLPSRLLRRIVHMEHQEADGDCHLQAVVL
HLARRSVCVHPQNRSLARWLERQGKRLQGTVP SLNLVLQKMYSHPPQQN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001199965



[View online »](#)

ORF Size:	360 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001199965.1</u> , <u>NP_001186894.1</u>
RefSeq Size:	425 bp
RefSeq ORF:	363 bp
Locus ID:	100039939
UniProt ID:	<u>Q9Z1X0</u>
Cytogenetics:	4
MW:	13.5 kDa