

## Product datasheet for **MR227942A1V**

### Mouse **Ccl27b (NM\_001199959) AAV Particle**

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Ccl27b (NM\_001199959) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Ccl27b  
**Synonyms:** Ctack; ENSMUSG00000066202; Scya27b  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR227942 representing NM\_001199959  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGGATCGCC

ATGATGGAGGGGCTCTCCCCGCCAGCAGCCTCCCGCTGTTACTGTTGCTTCTGAGCCCGGCTCCTGAAG  
 CAGCCTTGCTCTGCCCTCCAGCACTAGCTGCTGTACTCAGCTCTATAGACAGCCACTCCCAAGCAGGCT  
 GCTGAGGAGGATTGTCCACATGGAATGCAGGAGGCTGATGGGGACTGTCACCTCCAGGCTGTCGTGCTT  
 CACCTGGCTCGGCGCAGTGTCTGTGTTCCATCCCAGAACCCGAGCCTGGCTCGGTGGTTAGAACGCCAAG  
 GGAAAAGGCTCCAAGAACTGTACCCAGTTTAAATCTGGTACTACAAAAGAAAATGTACTCACACCCCA  
 ACAGCAAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR227942 representing NM\_001199959  
Red=Cloning site Green=Tags(s)

MMEGLSPASSLPLLLLLSPAPEAALPLPSSTSCCTQLYRQPLPSRLRRIVHMELQEADGDCHLQAVVL  
 HLARRSVCVHPQNRSLARWLERQGKRLQGTVP SLNLVLQKMYSHPPQQN

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_001199959



[View online »](#)

<b>ORF Size:</b>	360 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001199959.1</a></u> , <u><a href="#">NP_001186888.1</a></u>
<b>RefSeq Size:</b>	425 bp
<b>RefSeq ORF:</b>	363 bp
<b>Locus ID:</b>	100040048
<b>UniProt ID:</b>	<u><a href="#">Q9Z1X0</a></u>
<b>Cytogenetics:</b>	4
<b>MW:</b>	13.5 kDa