

Product datasheet for **MR227960A1V**

Mouse **Lyrm1 (NM_001285959) AAV Particle**

Product data:

Product Type: AAV Particles
Product Name: Mouse Lyrm1 (NM_001285959) AAV Particle
Tag: Myc-DDK
Symbol: Lyrm1
Synonyms: 1110065L10Rik; 2310004B22Rik; 4930404J24Rik
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227960 representing NM_001285959
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACAGCAGCAACCCGACAAGAAGTCCTCTGCCTCTACCGAAGCATTTCAGGCTTGCACGGAAATGGC
AGGCAGCATCAGGGCAGATGGAAGACACCATCGAAGAAAAACAGTACATACTAAAGGAAGCCAGAACACT
GTTCCAGAAAAACAAAATCTCACAGACCCAGAGCTGATTAACAATGCATAGTTGAATGCACGGCCCGG
ATTGAGATTGGACTGCATTACCAGATCCCTTATCCACGCCCTATCCATCTGCCTCCAATGGGTCTACCCG
CAGCAGAGGTGGGGACTTCAAACCCAGGAAAACTGAGAAAACCTTTCAAACCAAGTGTATCTGAAGTC
TCATGATGAAGTTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227960 representing NM_001285959
Red=Cloning site Green=Tags(s)

MTAATRQEVLCLYRSIFRLARKWQAASQMEDTIEEKQYILKEARTLFQKNKNTDPELIKQCI VECTAR
IEIGLHYQIPYPRPIHLPPMGLTARRRGLQTQEKLRLKSKPVYLKSHDEVS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001285959



[View online »](#)

ORF Size:	369 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001285959.1</u> , <u>NP_001272888.1</u>
RefSeq Size:	898 bp
RefSeq ORF:	369 bp
Locus ID:	73919
UniProt ID:	<u>Q9CQB7</u>
Cytogenetics:	7 F2
MW:	14.2 kDa