

Product datasheet for **MR227976A1V**

Mouse Rpl31 (NM_001252218) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Rpl31 (NM_001252218) AAV Particle
Tag: Myc-DDK
Symbol: Rpl31
Synonyms: Dts
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227976 representing NM_001252218
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCCCGCAAAGAAGGGTGGCGAGAAGAAGAAGGGCCGTTCTGCCATCAACGAGGTGGTGACCCGAG
AATACACCATCAACATTCACAAGCGCATCCATGGAGTGGGCTTCAAGAAGCGTGCTCCTCGGCACTCAA
AGAAATTCGGAAGTTTGCCATGAAGGAAATGGGACACCAGATGTGCGCATTGACACCAGGCTCAATAAA
GCCGTCTGGCCAAGGGAATAAGGAACGTTCCATATCGCATCCGAGTACGCTTGCCAGAAAACGTAATG
AGGATGAGGATCCCCAAACAAGCTCTACACGCTGGTAACTTACGTGCCTGTTACCACATTCAAAAATCT
ACAGACGGTCAACGTGGATGAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227976 representing NM_001252218
 Red=Cloning site Green=Tags(s)

MAPAKKGGEKKKGRSAINEVVTREYTIINHKRIHGVGFKRAPRALKEIRKFAMKEMGTPDVRIDTRLNK
AVWAKGIRNVPYRIRVRLSRKRNEDESPNKLYTLVTVYVPVTTFKNLQTVNVNEN

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001252218



[View online »](#)

ORF Size:	378 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001252218.1</u> , <u>NP_001239147.1</u>
RefSeq Size:	1261 bp
RefSeq ORF:	378 bp
Locus ID:	114641
UniProt ID:	<u>P62900</u>
Cytogenetics:	1 B
MW:	14.5 kDa