

Product datasheet for **MR227984A1V**

Mouse H2bc6 (NM_001290530) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse H2bc6 (NM_001290530) AAV Particle
Tag: Myc-DDK
Symbol: H2bc6
Synonyms: Gm11398; H2bc4; H2bc8; Hist1h2; Hist1h2be
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227984 representing NM_001290530
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGCCTGAGCCAGCCAAGTCCGCTCCCGCCCCGAAGAAGGGCTCCAAGAAGGCTGTCACCAAGGCCCAGA
 AGAAGGACGGCAAGAAGCGCAAGCGCAGCCGCAAGGAGAGCTACTCGGTGTACGTGTACAAGGTGCTGAA
 GCAAGTGCACCCCGACACCGGCATCTCCTCCAAGGCCATGGGCATCATGAACTCGTTCGTGAACGACATC
 TTCGAGCGCATCGCGGGCGAGGGCTCCCGCTGGCGCATTACAACAAGCGCTCGACCATCACATCCCGGG
 AGATCCAGACGGCCGTGCGCTGCTGCTGCCCGGGGAGCTGGCCAAGCACGCGGTGTCGGAGGGCACCAA
 GGCTGTCACCAAGTACACCAGCTCCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAAGTTTAA

Protein Sequence: >MR227984 representing NM_001290530
Red=Cloning site Green=Tags(s)

MPEPAKSAPAPKKGSKKAVTKAQKDKGKRRKRSRKESYSVYVYKVLKQVHPDTGISSKAMGIMNSFVNDI
 FERIAGEASRLAHYNKRSTITSREIQTAVRLLLPGELAKHAVSEGTKAVTKYTSSK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001290530



[View online »](#)

ORF Size:	378 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001290530.1</u> , <u>NP_001277459.1</u>
RefSeq Size:	445 bp
RefSeq ORF:	381 bp
Locus ID:	319179
UniProt ID:	<u>Q6ZWY9</u>
Cytogenetics:	13 A3.1
MW:	13.9 kDa