

Product datasheet for **MR227992A1V**

Mouse 3830406C13Rik (NM_001284385) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse 3830406C13Rik (NM_001284385) AAV Particle
Tag: Myc-DDK
Symbol: 3830406C13Rik
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227992 representing NM_001284385
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTTCCTTGTTTACTCAAGAAATTCACCTTTCTAAAAGACATGAAGAAATATTATCACAAAGATTAA
 TGTTACTTCAGAAAATGAAAACAAATTTGGAGATGAAAACACAGAAAGGCATCTCCTCCAAGCCAC
 TGAGACTGCTTCTAGAAGAACTCCGACTCTGAAGGACATAGATGCTGCAGAGGAGGCTTTTCAGACC
 AAGCTTATCCACACCCACAGCCTTCGATGCTTTCTCTTGAGACTCGTTACTGGGCATCAGTAGAAGAAC
 ATATCCCCAAGTGGGAAGTGTTCCTTTAGGAAGAGCACCATATCCAATTGGTGCTGAAAATCAAATGA
 GGTACCATTGTTCAAACCGAGGCAAAGTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227992 representing NM_001284385
 Red=Cloning site Green=Tags(s)

MASLFTQEIHLKRHEEILSQRLMLLQKMKNFQDENTERASLLQATETASRRNLRLLKIDAAEEAFQT
 KLIPHPQPSMLSLETRYWASVEEHIPKWELFLLGRAPYPIGAENQNEVPFVQTEAKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001284385
ORF Size: 384 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001284385.1, NP_001271314.1</u>
RefSeq Size:	2498 bp
RefSeq ORF:	384 bp
Locus ID:	218734
UniProt ID:	<u>Q8BGD0</u>
Cytogenetics:	14 A1
MW:	14.7 kDa