

Product datasheet for **MR228009**

Cds2 (NM_001291040) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cds2 (NM_001291040) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cds2
Synonyms: 5730450N06Rik; 5730460C18Rik; AI854580; D2Wsu127e
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228009 representing NM_001291040
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGACCGAACTACGGCAGAGGGTTGTCCGCGAGGACGCACCTCCGGAGGACAAGGAATCAGAGTCCGAAG
 CAAAGCTAGATGGAGAGACAGCGTCTGACAGCGAGAGCCGAGCGGAAACAGCTCCCCTACCAACCTCCGT
 AGATGACACCCCAGAGGTCCTTAACAGGGCCCTTTCCAACCTGTCCTCAAGGTGGAAGAACTGGTGGGTG
 AGAGGCATCCTGACTTTGGCCATGATCGCGTTTTTCTTCATTATCATTTACCTGGGACCAATGGTTTTGA
 TGATGATTGTTATGTGTGCCAGATTAAGTGTTCATGAAATAATCACTATTGGCTACAATGTATACCA
 CTCTACGACCTGCCCTGGTTCAGGACCTCAGCTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228009 representing NM_001291040
Red=Cloning site Green=Tags(s)


MTELQRVVRDAPPEDKESESEAKLDGETASDESRAETAPLPTSVDDTPEVLNRALSNLSSRWKNWWV
 RGILTAMIAFFFIILGPMVLMIMVCVQIKCFHEIITIGYNYHSYDLPWFRTLSW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM 001291040

ORF Size: 387 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001291040.1](#), [NP_001277969.1](#)

RefSeq Size: 1261 bp

RefSeq ORF:	390 bp
Locus ID:	110911
Cytogenetics:	2 64.15 cM
MW:	15.4 kDa
Gene Summary:	Provides CDP-diacylglycerol, an important precursor for the synthesis of phosphatidylinositol, phosphatidylglycerol, and cardiolipin.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228009