

## Product datasheet for **MR228027A1V**

### Mouse Anapc15 (NM\_001291352) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Anapc15 (NM\_001291352) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Anapc15  
**Synonyms:** 3200002M19Rik; 6330414C15Rik; APC15  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR228027 representing NM\_001291352  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGGATCGCC

ATGGAGCCAGGACTCAAACTCAGCTTAGGAGCCATGTCCACCTTGTCCCGTCACTCTTCCCTCGTGTGA  
 CTGAGACACTGTGGTTAATCTGGACCGACCCTGTGTGGAGGAGACAGAGCTGCAGCAGCAGGAGCAGCA  
 GCATCAGGCCTGGCTCCAAGCATCGCAGAGAAAGACAACAACCTGGTACCAATTGGCAAGCCGGCCTCA  
 GAGCACTACGATGATGAGGAAGAAGAGGATGATGAAGATGATGAGGACAGTGAAGAGGATTCCGAAGATG  
 ATGAGGACATGCAAGACATGGATGAGATGAATGACTATAATGAGTCACCTGATGATGGAGAGGTCAATGA  
 GGTGGACATGGAAGGCAACGAACAGGATCAGGACCAGTGGATGATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR228027 representing NM\_001291352  
Red=Cloning site Green=Tags(s)

MEPGLKLSLGAMSTLFPSTLFPRTETLWFNLDRPCVEETELQQQEQQHQAWLQSI AEKDNNLVPIGKPAS  
 EHYDDEEEDEDEDEDESEEDSEDEDMQDMEMNDYNESPDDGEVNEVDMEGNEQDQDQWMI

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_001291352



[View online »](#)

<b>ORF Size:</b>	399 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001291352.1</a></u> , <u><a href="#">NP_001278281.1</a></u>
<b>RefSeq Size:</b>	843 bp
<b>RefSeq ORF:</b>	399 bp
<b>Locus ID:</b>	75430
<b>UniProt ID:</b>	<u><a href="#">P60007</a></u>
<b>Cytogenetics:</b>	7 E2
<b>MW:</b>	15.4 kDa