

Product datasheet for **MR228034**

Serf2 (NM_001290837) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Serf2 (NM_001290837) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Serf2
Synonyms: C80083; m4F5rel; Msmac1l
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228034 representing NM_001290837
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGACCCGCGGTAACAGCGAGAGCTCGCCGCCAGAAGAACATGAAGAAGCAGAGCGACTCGGTTAAGG
GAAAGCGCCGAGATGATGGGCTTCTGCTGCCGCCGCAAGCAGAGTGCCCATCATCTCTACCCCCAGG
GACTCGGAGATCATGCAGCAGAAGCAGAAAAAGGCAACGAGAAGAAGGAGGAACCAAGTAGCCTTGTC
GCTTCGTGTCCAACCCTCTGCCCTCCGCTGTGTGCCTGGAGCCAGTCCCACCATGCTCGAGTTTCTTC
CTGTAGTGCTCACAGGTCCAGCACCGATGGCATTCCCTTTGCCCTGAGTCTGCAGCGGGTTCCTTTTGT
GCTTCCTTCCCCTCAGGTAGCCTCTCTCCCTCTGGGCCACTCCTGGGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228034 representing NM_001290837
Red=Cloning site Green=Tags(s)

MTRGNQRELARQKNMKQSDSVKGKRRDDGLSAAARKQSAPSSLPPGTRRSCSRSRKRQTRRRRNPSLV
ASCPTLLPSACVPGASPTMLEFLPVVLTGPSTDGIPFALSLQRPVFLPSPQVASLPLGHSWG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001290837

ORF Size: 399 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001290837.1](#), [NP_001277766.1](#)

RefSeq Size: 3095 bp

RefSeq ORF:	402 bp
Locus ID:	378702
Cytogenetics:	2 E5
MW:	14.8 kDa
Gene Summary:	Positive regulator of amyloid protein aggregation and proteotoxicity (By similarity). Induces conformational changes in amyloid proteins, such as HTT, driving them into compact formations preceding the formation of aggregates (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228034