

## Product datasheet for **MR228045**

### Cck (NM\_001284508) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cck (NM\_001284508) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Cck  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR228045 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGAAGAGCGGCGTATGTCTGTGCGTGGTATGGCAGTCCTAGCTGCTGGCGCCCTGGCGCAGCCGGTAG  
 TCCCTGCAGAAGCTACGGACCCCGTGGAGCAGCGGGCGCAAGAGGCGCCCGAAGGCAGCTGCGGGCTGT  
 GCTCCGGACGGACGGCGAGCCCCGAGCGCGCTGGGCGCACTGCTAGCGCGATACATCCAGCAGGTCCGC  
 AAAGTGGCATGGATGGTACCTCTGGTTGGGTGTTAACCTGGACCAGTAGAGCTGGATTTAAACACAGGA  
 GGTGGGCATCCTTTCTGTGGAGCTCACGAACCAATTTTCTGCCCGCATTGTAACAGCCCATGGCTTG  
 TCGACCTGTCTGCATTTGGCTTGACTGCTCCTTCTGGCCGCATGTCCGTTCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR228045 protein sequence  
 Red=Cloning site Green=Tags(s)

MKSGVCLCVMAVLAAGALAQPVVPAEATDPVEQRAQEAPRRQLRAVLRDGEPRARLGALLARYIQQVR  
 KVAWMVTSGWVLTWTSRAGLKHRRWASFLWSSRTQFFLPAFEQPMACRPVCIWLDCSFVPHVRS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001284508

**ORF Size:** 405 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001284508.2](#), [NP\\_001271437.1](#)


**RefSeq Size:** 854 bp

**RefSeq ORF:** 405 bp

**Locus ID:** 12424  
**UniProt ID:** [P09240](#)  
**Cytogenetics:** 9 F4  
**MW:** 15.2 kDa

**Gene Summary:** This gene encodes a member of the gastrin/cholecystokinin family of proteins. The encoded preproprotein is proteolytically processed to generate multiple protein products, including the peptide hormones cholecystokinin-8, -12, -33, and others. The encoded peptides have been shown to regulate gastric acid secretion and food intake. A sulfated form of cholecystokinin-8 may modulate neuronal activity in the brain. Homozygous knockout mice for this gene exhibit impaired insulin secretion, enhanced insulin sensitivity, and resistance to obesity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2015]

**Product images:**



Circular map for MR228045