

Product datasheet for MR228092

Csmd2os (NM_029137) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Csmd2os (NM_029137) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Csmd2os
Synonyms: 4930465A12Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228092 representing NM_029137
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGAGAAGAAGCAGCAAGGAACCTCATGGGTTGGATGAGGCCAGCTCTGCTCCATGGGAGTCAAGG
 TGATGGGGAGAAAGTTGAGTGATGACCCAAAGACGTTTGCAATTCTAGCAGTGTCATGGCTCAGCTGC
 AGCCCAGGAAGTCCCAGGTCCTGGGTACATAGCGAGGGAGGAGCATTGATCTGCTAGAAGCTGTGCCA
 AGGTCCGTGGGAACCGCAACCCCTTCCCTCGACATCTTGCCAGCAGCACTGGATTAATCTCTGAGAACA
 AGCACTCTGGGTTCCAGAGGCAAAACAGGCTTTTGATTCTATGAGCTAAGACTCCTCTGCAGGAGCGT
 GGTCTTGTGCCTTACTGTCATTGCTTCTTACAAGGATCACATAAGCCAAAGATCCCAGGAACCTCCTGTG
 GGCAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228092 representing NM_029137
 Red=Cloning site Green=Tags(s)

MEKKQQGTPHGLDEAQLCSMGVKVMGRKLSDDPKDVCISSSVHGSAAAQEVPRSWVHSDGRSIDLLEAVP
 RSVGTATSPRHLASSTGLISENKHSGFPEAKQAFDFYELRLCRSVVLCLTVIASYKDHISQRSQEPVP
 GK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_029137

ORF Size: 426 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_029137.1](#), [NP_083413.1](#)


RefSeq Size: 643 bp

RefSeq ORF: 429 bp

Locus ID: 74989

Cytogenetics: 4 D2.2
MW: 15.8 kDa

Product images:



Circular map for MR228092