

## Product datasheet for **MR228094**

### Elk3 (NM\_001282967) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Elk3 (NM\_001282967) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Elk3  
**Synonyms:** D430049E23Rik; Erp; Etrp; Net; Sap-2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR228094 representing NM\_001282967  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGAGAGTGAATCACGCTGTGGCAGTTCCTCTTGCACTTGCTGCTGGACCAGAAACATGAGCACCTCA  
TCTGCTGGACATCGAACGATGGCGAGTTCAAGCTCCTCAAGCGGAAGAAGTGGCCAAGCTGTGGGGCT  
CCGCAAGAACAAGACCAACATGAACTACGACAAGCTGAGCAGAGCGCTGAGATACTATTACGACAAGACA  
CCAAGTGGACTGTTTCTGGCCTCGAGTCCGCTGCTGCCAGCATACACTTCTGGAGCAGCCTTAGTCCTG  
TGCCCCACTGAGTCTGCCAGGCTGCAAGGGCCGAACACACTTTTCCAGTTCGCCACACTGCTCAACGG  
TCACATGCCGGTGCCGCTCCCCAGTCTGGACAGAGCTCCATCCCCAGTTCGCTGTCCCCAGCTCTCAG  
AAATCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR228094 representing NM\_001282967  
Red=Cloning site Green=Tags(s)

MESAITLWQFLLHLLLDQKHEHLICWTSNDGEFKLLKAEVAKLWGLRKNKTNMNYDKLSRALRYYYDKT  
PSGLFLASSPLLPSIHFWSLSPVAPLSPARLQGPNTLFQFPTLLNGHMPVPLPSLDRAPSPVLLSPSSQ  
KS


**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001282967

**ORF Size:** 426 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001282967.1, NP_001269896.1</u>
<b>RefSeq Size:</b>	3433 bp
<b>RefSeq ORF:</b>	429 bp
<b>Locus ID:</b>	13713
<b>Cytogenetics:</b>	10 48.04 cM
<b>MW:</b>	16.5 kDa
<b>Gene Summary:</b>	May be a negative regulator of transcription, but can activate transcription when coexpressed with Ras, Src or Mos. Forms a ternary complex with the serum response factor and the ETS and SRF motifs of the Fos serum response element.[UniProtKB/Swiss-Prot Function]