

Product datasheet for MR228097

Prr3 (NM_001282013) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prr3 (NM_001282013) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Prr3
Synonyms: 4930540G07Rik; Cat56
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228097 representing NM_001282013
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTAACGGGAAGCCTGGTGATCCCAAGTCAGCTTTCCACAGGGTCTCCAGGATCAAGGGACGGA
TGATCCCCACTGCTGAGTCTCCGCCTCTCCCGGGTAGAGGCTACATTCGGGAGGCCTCGGTCC
TCGATCTAGCCCGTATGGTCGTGGTTGGTGGGGTATCAATGCCGAACCTCCTTTCTGGCCAGGCCAT
GGGGTCCCTCTAGGGAAAGCTTTTACAAAGAGGCGAGAAACCCACGAAGGCTCAGAAGCTGGTCTCTTG
TCAAGAATACCTACCCACCAAGGACAGCCCCAGATGATGGAAGACAAGTCTGACCGCCCTGTCTGCCG
ACACTTTTCTAAGAAGGGCCACTGTCGATATGAAGACCATTGTGCCTTCTACCATCCAGGCGTCAACGGA
CCTCCTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228097 representing NM_001282013
Red=Cloning site Green=Tags(s)
MANGKPGDPKSAFHRGPPGSRGRMIPLLSLPPPPRGRGYIRGGLGPRSSPYGRGWGINAEPFPGPGH
GGPSRESFYKEARNPRRLRSWSLVKNTYPPKDSPOMMEDKSDRPVCRHF SKKGHCERYEDHCAFYPGVNG
PPL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001282013

ORF Size: 429 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001282013.2</u>
RefSeq Size:	1947 bp
RefSeq ORF:	432 bp
Locus ID:	75210
UniProt ID:	<u>Q811B5</u>
Cytogenetics:	17 B1
MW:	16.3 kDa