

Product datasheet for **MR228140**

1700015F17Rik (NM_001200025) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700015F17Rik (NM_001200025) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: 1700015F17Rik
Synonyms: Gm1686
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228140 representing NM_001200025
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTGCTGGCTGAGGGCATGGAGCCACATCCTCCTGCCGTTTTCTCTCTGTAGCTCTCATCCAGCTGA
TTTTCAACTTATCTGACCATAAGACTAAAACCCACCATCACTGGAAGCGGAAAATCGAAGCAAAATACGT
TGAAAAGGAATGTGCTCAAAAGGAAACATGTCCAGTGTGTAAGGATAAGAGATGATTTGGTGCCGA
GAAGAAAAGGTTTGCAAAAGTACTGTTTTCCATATTCTGATTGTAATTAATTCTATATTTGGGCAA
ACTGTAATGTTGACCTGTTTGAATCGTGATGCTGATACTGATTGTAATATTAGCACTAGCATTCTCTG
GTATTGCCTCGCCTATTACTTCTACATGCAACAGCATATGGCCCTTTATGCAAGACATGGACAAGTTCCC
GTCTACAACCTGGGATGCTCCAGGTGACTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228140 representing NM_001200025
Red=Cloning site Green=Tags(s)


MCWLRAWSHILLPVFLSVALIQLIFNLSDHKTKTHHHWKRKIEAKYVEKECAQKETCPVCTRDKRCIWC
 EEKVCKKYCFPYSDCKFNSIFWANCNVDLFGIVMLILIVILALFLWYCLAYFYMQQHMALYARHGQVP
 VYNWDAPGDW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001200025

ORF Size: 450 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001200025.1](#), [NP_001186954.1](#)

RefSeq Size: 710 bp

RefSeq ORF: 453 bp


Locus ID: 381716

UniProt ID: D3YUK8

Cytogenetics: 5 A1

MW: 18.5 kDa

Product images:



Circular map for MR228140