

Product datasheet for **MR228144**

Fkbp1a (NM_001302078) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fkbp1a (NM_001302078) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fkbp1a
Synonyms: Fkb; Fkbp; Fkbp1; FKBP12
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228144 representing NM_001302078
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACTTACAGGGAAGGGAAGACAAGACAGCCCGCTCACATCCTCCTCTGTGCCCCCTTACTTCTATG
ACTTCTTGAGTAGAGAAGCGGGCAGGAGATGCCAGAAGTTAAGGACTAGGACCTTTAAAAAGCTGCCAAT
GGAAGGAATATTGAAATACATCTTGGCAAGAACAGACTCTTCTCATGGACCAAATCGCAGTAGCTGTTGT
GGGATGCTTGAAGATGGAAAGAAATTTGATTCTCTCGGGACAGAAACAAGCCTTTTAAGTTTACTACTAG
GCAAGCAGGAGGTGATCCGAGGCTGGGAGGAAGGGGTAGCCAGATGAGTGTGGGTCAGAGAGCCAACT
GATAATCTCCTCAGACTATGCCTATGGAGCCACCGGGCACCCAGGCATCATCCCACCACATGCCACTCTT
GTTTTTGTGTTGGAGCTTCTAAACTGGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR228144 representing NM_001302078
Red=Cloning site Green=Tags(s)
MTYREGKTRQPAHILLSVPPYFYDFLSREAGRRWQKLRTTRTFKKLPMEGILKYILARTDSSHGPNRSSCC
GMLEDGKKFDSSRDNRNPKFKFTLGKQEVIRGWEEGVAQMSVQRAKLIISDYAYGATGHPGIIPPHATL
VFDVELLKLE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-Mlul



Cloning Scheme:

Plasmid Map:


ACCN: NM_001302078

ORF Size: 450 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001302078.1, NP_001289007.1</u>
RefSeq Size:	1902 bp
RefSeq ORF:	453 bp
Locus ID:	14225
UniProt ID:	<u>P26883</u>
Cytogenetics:	2 G3
MW:	17.6 kDa
Gene Summary:	This gene is a member of the immunophilin family. The encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin, and is associated with immunoregulation, protein folding, receptor signaling, protein trafficking and T-cell activation. It may modulate the calcium release activity of the ryanodine receptor Ryr1. It also interacts with the type I TGF-beta receptor. Disruption of this gene in mouse causes severe ventricular defects. Pseudogenes of this gene have been defined on chromosomes 4, 10, 14, and 16. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2014]