

Product datasheet for **MR228150**

Fam173a (NM_001285984) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam173a (NM_001285984) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fam173a
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228150 representing NM_001285984
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGGATCTAGGCTCTGGTGATGGCAGGATCGTGCTGGCGGCCACCAATGTGGTCTCCGTCCAGCTA
 TGGGATATGAGCTGAACCCGTGGCTTGTGGGACTTGCACGGCTGCATGCCTGGAGAGCCGGTTGTTCTGC
 GAGTGTCTGCTACCATCGAAAAGATCTCTGGAAGGTGAGCCTCAGGGACTGTCAACAGTGTCTGTGTTC
 CTGGCTCCCAAGTGTCTCAACTGCTGGAGGACAACTGCAAGCAGAGCTCCCTGTGGGGGCCCGTGTGG
 TATCTGGACGATTCCCCCTTCCACCTGGCAGCCTGTAGCTGTGGTGGGCGAGGGCACTGACAGAGTTTG
 GGCTATGATGTCCATGGTTCTGGGCCAACTGTGTCTTCTGTGGGTGCCCATCAAGGCCATCCCAGAG
 TCCAGTTCTACCTTAGTCCCTAGAGCCCCTGTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228150 representing NM_001285984
 Red=Cloning site Green=Tags(s)

MYDLGSGDGRIVLAAHQCLRPAMGYELNPWLVLARLHAWRAGCSASVCYHRKDLWKVSLRDCHNVSVF
 LAPSVLQLLEDKLAELPVGARVVSGRFLPTWQPVAVVGEGETDRVWAYDVHSGSPTVSSCGVPIKAIPV
 SSSTLVPRAPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001285984

ORF Size: 453 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001285984.1](#), [NP_001272913.1](#)

RefSeq Size: 633 bp

RefSeq ORF: 456 bp

Locus ID:	214917
Cytogenetics:	17 A3.3
MW:	16.7 kDa
Gene Summary:	S-adenosyl-L-methionin-dependent protein-lysine N-methyltransferase.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228150