

Product datasheet for MR228177

Ccdc50 (NM_001289436) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccdc50 (NM_001289436) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ccdc50
Synonyms: C3orf6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228177 representing NM_001289436
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAAACCAAGAGGAATAAAAGAAGCTGTCTCTACTCCAGCACGAGCAAGCCACAGGGACCAGGAGTGGT
ATGATGCTGAAATTGCCAGGAAATTGCAAGAAGAAGAACTTTTGGCTACTCATGTGGACATGAGAGCAGC
TCAGGTTGCCAGGATGAGGAAATTGCTCGACTTCTAATGGCTGAAGAAAAAAGCTTACAAGAAAGCC
AAAGAGCGAGAAAAGTCATCTTTGGACAAAAGGAAACATGACCCTGAATGCAAGTTAAAAGCAAAGTCAG
CCCCTCAAAGTCAAAGAGGGTGATGAAGCACACCGCTCCAAGATTGACAGGCCATCAAGACCACCACC
ACCTACCATGATGGGCCTTGAGGACACAGATCCCACCCATTTTACAAACCAGCACAGTACAACATGGCAT
CTTCAAAGTCAGAGTCCTCACAGAAAGGCTTCCATAACAAGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228177 representing NM_001289436
Red=Cloning site Green=Tags(s)

MKPRGIKEAVSTPARASHRDQEWYDAE IARKLQEEELLATHVDMRAAQVAQDEE IARLLMAEKKAYKKA
KEREKSSLDKRKHDPECKLKAKSAHSKSKEGDEAHRSKIDRPSRPPPTMMGLEDTDPHTFNQHSSTWH
LPKSESSQKGFHNKQ

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-Mlul




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001289436

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001289436.1, NP_001276365.1</u> |
| RefSeq Size: | 7063 bp |
| RefSeq ORF: | 468 bp |
| Locus ID: | 67501 |
| UniProt ID: | <u>Q810U5</u> |
| Cytogenetics: | 16 18.98 cM |
| MW: | 18.3 kDa |
| Gene Summary: | Involved in EGFR signaling.[UniProtKB/Swiss-Prot Function] |