

Product datasheet for MR228186

Hes2 (NM_001301805) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hes2 (NM_001301805) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Hes2
Synonyms: HES-; HES-2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228186 representing NM_001301805
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGGCTGCCCAGAAGAGTTGAGGACGCGCGGAGCTGCGCAAGAACCTAAAGCCGCTGCTGGAGAAGC
 GCCGGCGCGCGGATCAACGAGAGCCTAAGCCAGCTGAAGGGTCTCGTATTGCCGCTGCTGGGCGCGGA
 GACTTCCCCTCTTCGAAGTTGGAGAAGGCAGACATCCTGGAAATGACTGTGCGCTTCTACAGGAGCAG
 CCTGCGACCCTGTACTCCTCTGCAGCACCTGGCCTTTGAATAGCTACCTCGAGGGTTACCGAGCCTGTC
 TGGCGGACTGGCCCGCTACTGCCCGCTGCAGCGTCTGGAGCCCGCAGTGAGCGCGCGCCTACTAGA
 GCACCTGCGACAGAGGACTGTCAGTGATGACTCGCCCTACTGACGCTACCACCCGCTCCAGCTCCAGCT
 CCTTCTCCGCGGTGCCTCCTCCAGGCAGTTCTGGCCTTTGGAGGCCATGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228186 representing NM_001301805
 Red=Cloning site Green=Tags(s)

MRLPRRVEDAAELRKNLKP LLEKRRRARINESLSQLKGLVPLPGAETSRSSKLEKADILEMTVRFLQEQ
 PATLYSSAAPGPLNSYLEGYRACLARLARVLPACSVLEPAVSARLLEHLRQRTVSDDSPSLTLPPAPAPA
 PSPPVPPPGSSGLWRPW

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online >](#)

Cloning Scheme:


ACCN: NM_001301805

ORF Size: 471 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001301805.1](#), [NP_001288734.1](#)

RefSeq Size: 2915 bp

RefSeq ORF: 474 bp

Locus ID: 15206


UniProt ID: [O54792](#)

Cytogenetics: 4 82.92 cM

MW: 17.3 kDa

Gene Summary: The protein encoded by this gene belongs to the mammalian Hes gene family, the mammalian homologues of *Drosophila* hairy and Enhancer of split. Hes 2 is a basic helix-loop-helix transcriptional repressor and is an effector of the Notch signaling pathway. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2014]

Product images:



Circular map for MR228186