

Product datasheet for **MR228199**

Emp1 (NM_001288627) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Emp1
Synonyms:	I-8-09; TMP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR228199 representing NM_001288627
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGCTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGTTGGTGTCTACTGGCTGGTCTCTTTGTGGTCCACATTGCCACTGCCATTATGCTGTTGTCTCCACCA
TTGCCAACGTCTGGATGGTTGCAGATTACGAAATGCATCTGTAGGGCTTTGGAAGAACTGCACTGGTGG
TAATGCGACGGCTCCCTGTCTACGGCAATGAAGATGCTATCAAGGCAGTGCAAGCCTTCATGATCCTC
TCCATCATCTTCTCCATCATCTCCCTCGTGGTCTTCGTGTTCCAGCTCTTCACTATGGAGAAGGGAACC
GGTCTTCTCTCGGGGTCCACCATGCTGGTGTGCTGGCTGTGTATCCTGGTTGGAGTGCAATCTACAC
TCATCATTACGCCACAGCGAAGGGAACCTCAACTCCAGCAGCCACCAAGGCTATTGTTTCATCCTGACC
TGGATCTGCTTCTGTTTCAGCTTCATCATCGGCATACTCTATATGGTCCTGAGGAAGAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228199 representing NM_001288627
Red=Cloning site Green=Tags(s)

MLVLLAGLFVVHIATAIMLFVSTIANVWMVADYANASVGLWKNCTGGNCDGSLSYGNEDAIAKAVQAFMIL
SIFSIISLVVFVQLFTMEKGNRFFLSGSTMLVCWLCILVGVSIYTHHYAHSEGNFNSSSHQGYCFILT
WICFCFSFIIGILYMLRKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001288627

ORF Size: 480 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001288627.1](#), [NP_001275556.1](#)

RefSeq Size:	2826 bp
RefSeq ORF:	483 bp
Locus ID:	13730
UniProt ID:	P47801
Cytogenetics:	6 66.25 cM
MW:	17.8 kDa

Product images:



Circular map for MR228199