

Product datasheet for **MR228215**

Ppil3 (NM_001285826) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Ppil3
Synonyms:	2310076N22Rik; 2510026K04Rik; Cyp10l
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR228215 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCTGTGACACTGCATACAGATGTGGGAGATATCAAGATAGAAGTCTTCTGTGAGAGAACACCCAAAA
CATGTGAGAATTTCTTGCTCTTTGTGCCAGTAATTACTACAACGGCTGTGTGTTTCATAGAAATATCAA
GGGCTTCATGGTTCAAACAGGAGATCCAACAGGTACTGGAAGAGGAGGTAGCAGCATCTGGGCCAAAAAG
TTTGAGGATGAATACAGTGAATATCTGAAGCACAATGTTGAGGTGTTGTATCTATGGCTAATAATGGCC
CAAATACCAATGGATCTCAGTTCTTCATCACCTATGGCAAGCAGCCACACTTGGACATGAAATATACAGT
GTTTGAAAGGTAATAGATGGTCTGGAGACTTTGGATGAGTTGGAGAAAGTTACCAAGTAAATGAGAAGACA
TACAGACCTCTTAATGATGTACACATTAAGGACATAACTATTCATGCCAACCCATTTGCTCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228215 protein sequence
Red=Cloning site Green=Tags(s)

MSVTLHTDVGDIKIEVFCERTPKTCENFLALCASNYNGCVFHRNIKGFVMVQTGDPTGTGRGSSIWAKK
FEDEYSEYLLKHNVRGVVSMANNPNTNGSQFFITYGKQPHLDMKYTVFGKVIDGLETLDLEKLVPNEKT
YRPLNDVHIKDITIHPFAQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001285826

ORF Size: 483 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001285826.1](#), [NP_001272755.1](#)

RefSeq Size: 1421 bp

RefSeq ORF: 486 bp

Locus ID: 70225


UniProt ID: [Q9D6L8](#)

Cytogenetics: 1 C1.3

MW: 18.1 kDa

Gene Summary: PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. May be involved in pre-mRNA splicing (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228215