

Product datasheet for MR228252

Zfp740 (NM_001289691) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Zfp740 (NM_001289691) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Zfp740
Synonyms: 1110034O07Rik; AW548228; Znf740
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228252 representing NM_001289691
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGATGCTGAGCCAGATTGCCAGCAAGCAGGCTGAAAACGGCGAGCGGGCAGGTAGCCCTGATGTGCTGA
 GGTGCTCCAGTCAGGGCCACCGAAAAGACAGTGATAAATCCCGCAACCGCAAAGAGGATGACAGCTTGCC
 TGAGGCCTCTCATTCAAAAAAGACTGTTAAAAAGGTGGTGTAGTGAACAAAATGGCTCTTTCAAGTA
 AAGATTCCTCAAAAATTTTATTTGTGAACACTGCTTTGGAGCCTTCGGAGCAGTTACCACCTCAAGAGGC
 ACGTCCTAATCCCACTGGTGAAGAACATTTGAGTGTGATGTCTGTGATATGCGTTTCATCCAGAAGTA
 CCATCTCGAACGCCACAAGCGGGTACACAGTGGCGAAAAGCCTTACCAGTGTGAACGATGTCATCAGTGT
 TTTTCTCGGACAGACCGATTACTCAGACACAAACGGATGTGCCAAGGATGCCAGTCCAAGACTTCTGAAG
 GGCAGTTTTCTCTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228252 representing NM_001289691
 Red=Cloning site Green=Tags(s)

MMLSQIASKQAENGERAGSPDVLRCSSQGHRKDSKSRNRKEDDSLAEASHSKTVKKVVVVEQNGSFQV
 KIPKNFICEHCFGAFRSSYHLKRHLIHTGEKPFECVCDMRFIQKYHLERHKRVHSGEKPYPQERCHQC
 FSRTDRLLRHKRMCQGCQSKTSEQQFSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001289691

ORF Size: 504 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001289691.1](#), [NP_001276620.1](#)


RefSeq Size: 4249 bp

RefSeq ORF: 507 bp

Locus ID: 68744

UniProt ID: Q6NZQ6
 Cytogenetics: 15 F2-F3
 MW: 19.9 kDa
 Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228252