

Product datasheet for **MR228258**

Ssxb6 (NM_001205108) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ssxb6 (NM_001205108) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ssxb6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228258 representing NM_001205108
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAACAGTGAGTTCCTGTGAAAAGATCCCCATGGAAGTGCTTTATGAACAAAGAATATCTGCAAGG
CATTCCAGGATATTTCTACATACTTCTCTGATGAAGAGTGGTAAAGCTGACTCAATGGCAGAAAAGTGC
CTATGTGTACATGAAAAGAACTACATCAGAATGACTGACCTAGGTGTCACTGTGAACCAACCAGTTTC
ATGCGTGGCAAGGAGCAGGCCAAGCAATCCCTGGTGAAGGCATTGAAATTCATGACAGTGATGATGAAT
GCTTTGAAGGATCTTTTGGTGTGACACCGATAAAACGAATGAAGCTGACATCAGTGACAATAAGTATTCA
TAATGTAGAAGGAAGCCATGACTCTGGAGAGAATGATTCAAACCTGGCAGAACTGGTGGCATTACGGTC
AATGTTTGGAGCCACCGCTGCGAGAGAGGAAGTATCGAGTGATATATTCAGAGATTAGTGATCCTGAAG
AGGAAGAAGATGATGATTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228258 representing NM_001205108
Red=Cloning site **Green**=Tags(s)


METVSSCEKIPMEVL YEPKNICKAFQDISTYFSDEEWVKLTQWQKSAYVYMKRNYIRMTDLGVTNVQPVF
MRGKEQAKQSLVEGIEIHDSDECFEGSFGVTPIKRMKLSVTISIHNVESGSHDSGENDSNLAETGGIQV
NVWSHRLRERKYRVIYSEISDPEEEEDDDY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001205108

ORF Size: 510 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001205108.1](#), [NP_001192037.1](#)

RefSeq Size: 522 bp

RefSeq ORF: 513 bp

Locus ID: 668976
 Cytogenetics: X A1.1
 MW: 20.2 kDa

Product images:



Circular map for MR228258