

Product datasheet for **MR228301**

Styx11 (NM_001289556) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Styx11 (NM_001289556) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Styx11
Synonyms: 1700011C14Rik; Dusp24; Mkstyx
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228301 representing NM_001289556
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTACCACTTCTCCGGACTCAGAAGATCATCTGGATGCCCCAGGAAGTGGATGCCTTCCAGCCATACC
CAGTTGAGATACTCCAGGCAGAGTCTACCTGGGCAAAATCAGTCAAGCCTGCAATGCTAAGATGCACAA
GGACTTGAAAATTAAGCACATGCAATATTTCCATGGAGACAACACCCTATTTTATAGGCAATGCAGAC
AAGCTCTTACATCAAGCTTGAGGATTCTCCAGACACCTTGCTTTTGGACTATTTCCGCCACATCTGTC
ACTTTATAGAAGTTCACCTCAGGCTCCACTCTGTCGTTCTGGTCTTTTCCACCCGGGCATCAGCCGCAG
CGTTGCTGCTGTGGTGGCCTTACTCATGCATTATAATGAGGAGACGCTGAAGAGGTCTGGGCCAAAGTG
AAGAAATGCAAAACCAACATGCGTCCGAATCGGGCCTTCGTGGCTCAGCTGCTGGAGTGGGAGAAAGTCC
TCCATGGGGAATATCTCACAGACATCTCCGAGCCTATCTAC

ACGCGTACGCGGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228301 representing NM_001289556
Red=Cloning site Green=Tags(s)

MYHFFRTQKIIWMPQELDAFQYPVEILPGRVYLGKISQACNAKMHKDLKIKAHVNISETTPYFIGNAD
KLLHIKLEDSPDTLLFDYFRHICHFIELHLRLHSVVLVSTRGISRSVAVVALLMHYNEETLKRSWAQV
KKCKTNMRPNRAFVAQLLEWEKVLHGEYLTDISEPIY


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online >](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001289556

ORF Size: 531 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001289556.1</u> , <u>NP_001276485.1</u> |
| RefSeq Size: | 1511 bp |
| RefSeq ORF: | 534 bp |
| Locus ID: | 76571 |
| Cytogenetics: | 5 G2 |
| MW: | 21.2 kDa |