

Product datasheet for MR228321

Zfp740 (NM_001289695) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Zfp740 (NM_001289695) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Zfp740
Synonyms: 1110034O07Rik; AW548228; Znf740
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228321 representing NM_001289695
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGATGCTGAGCCAGATTGCCAGCAAGCAGGCTGAAAACGGCGAGCGGGCAGGTAGCCCTGATGTGCTGA
 GGTGCTCCAGTCAGATGGACTGTAAGCCTAGATTTGATTTGTCTTCAAAGGGCCACCGAAAAGACAGTGA
 TAAATCCCGCAACCGCAAAGAGGATGACAGCTTGGCTGAGGCCTCTCATTCAAAAAGACTGTAAAAAG
 GTGGTGGTAGTGAACAAAATGGCTCTTTTCAAGTAAAGATCCCAAAAATTTTATTTGTGAACACTGCT
 TTGGAGCCTTTCCGAGCAGTTACCACCTCAAGAGGCACGTCTTAATCCCACTGGTGAGAAACCATTGGA
 GTGTGATGTCTGTGATATGCGTTTCATCCAGAAGTACCATCTCGAACGCCACAAGCGGGTACACAGTGGC
 GAAAAGCCTTACCAGTGTGAACGATGTCATCAGTGTTTTTCTCGGACAGACCGATTACTCAGACACAAAAC
 GGATGTGCCAAGGATGCCAGTCCAAGACTTCTGAAGGGCAGTTTTCTCTA

ACGCGTACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228321 representing NM_001289695
 Red=Cloning site Green=Tags(s)


MMLSQIASKQAENGERAGSPDVLRCSSQMDCKPRFDLSSKGHRKDSKSRNRKEDDSLAEASHSKKTVKK
 VVVVEQNGSFQVKIPKNFICEHCFGAFRSSYHLKRHLIHTGEKPFECDCVDMRFIQKYHLERHKRVHSG
 EKPYQCERCHQCFSRDRLLRHKRMCQGCQSKTSEGQFSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:

Plasmid Map:


ACCN: NM_001289695

ORF Size: 540 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289695.1</u> , <u>NP_001276624.1</u>
RefSeq Size:	3924 bp
RefSeq ORF:	543 bp
Locus ID:	68744
UniProt ID:	<u>Q6NZQ6</u>
Cytogenetics:	15 F2-F3
MW:	21.3 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]