

Product datasheet for MR228322

Mup9 (NM_001281979) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mup9 (NM_001281979) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Mup9
Synonyms: Gm14076; Mup2; Mup6; OTTMUSG00000015595
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228322 representing NM_001281979
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGAAGATGCTGCTGCTGCTGTGTTGGGACTGACCCTAGTCTGTGTCCATGCAGAAGAAGCTAGTTCTA
CGGGAAGGAACCTTAATGTAGAAAAGATTAATGGGGAATGGCATACTATTATCCTGGCCTCTGACAAAAG
AGAAAAGATAGAAGATAATGGCAACTTTAGACTTTTTCTGGAGCAATCCATGTCTTGAGAATTCCTTA
GTTCTTAAATTCATACTGTAAGAGATGAAGAGTGCTCCGAATTATCTATGTTGCTGACAAAACAGAAA
AGGCTGGTGAATATTCTGTGACGTATGATGGATTCAATACATTTACTATACCTAAGACAGACTATGATAA
CTTTCTTATGGCTCATCTCATTAAACGAAAAGGATGGGGAACCTTCCAGCTGATGGGGCTCTATGGCCGA
GAACCAGATTTGAGTTCAGACATCAAGGAAAGGTTTGCAAACTATGTGAGGAGCATGGAATCCTTAGAG
AAAATATCATTGACCTATCCAATGCCAATCGCTGCCTCCAGGCCGAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228322 representing NM_001281979
Red=Cloning site Green=Tags(s)

MKMLLLCLGLTLVCHAEASSTGRNFNVEKINGEWHTIILASDKREKIEDNGNFRFLFLEQIHVLENSL
 VLKFHTVRDEECSEL SMVADKTEKAGEYSVTYDGFNTFTIPKTDYDNFLMAHL INEKDGETFQLMGLYGR
 EPDLSSDIKERFAQLCEEHGILRENIIDLSNANRCLQARE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM 001281979

ORF Size: 540 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001281979.1](#), [NP_001268908.1](#)

RefSeq Size: 882 bp

RefSeq ORF: 543 bp


Locus ID: 100038948

UniProt ID: P02762

Cytogenetics: 4

MW: 21.1 kDa

Product images:



Circular map for MR228322