

Product datasheet for **MR228364**

Vti1a (NM_001293686) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vti1a (NM_001293686) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Vti1a
Synonyms: 1110014F16Rik; 1110018K19Rik; 4921537J05Rik; MVti1; MVti1a; Vti1; Vti1-rp2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228364 representing NM_001293686
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTTCGACTTCGAAGGGTACGAGCAGGACTTCGCAGTGTCTACTGCTGAGATCACCAGCAAGATTG
CAAGGGTCCCTCGCTCCCTCCCGATGAGAAGAAGCAAATGGTTGCCAATGTGGAGAAACAGCTTGAAGA
GGCCAGAGAACTGCTTGAACAGATGGATTTGGAAGTTCGAGAAATCCACCCCAAAGTCGAGGAATGTAT
AGCAACAGGATGAGAAGCTACAAGCAAGAAATGGGGAAACTGGAAACGGATTTTAAAAGGTCACGGATCG
CCTACAGTGACGAAGTACGGAATGAGCTCTAGGGGATGCTGGGAATTCCTCCGAGAACCAGAGGGCACA
TCTGCTGGATAACACGGAGAGGCTGAAAGGTCGTCTCGGAGACTAGAGGCTGGGTACCAGATAGCAGTG
GAGACTGAGCAAATCGGCCAAGAGATGTTGGAAAACCTGAGTCATGACAGAGAAAAGATACAGCGGGCAC
GTGACCGACTCCGGGACGCAGATGCAAACCTGGGAAGAGCTCTCGGATTCTGACGGGGATGCTCGGAAG
G

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228364 representing NM_001293686
Red=Cloning site Green=Tags(s)



MSSDFEGYEQDFAVLTAEITSKIARVPRLPPEDEKKQMVANVEKQLEEARELLEQMDLEVREIPPQSRGMY
SNRMRSYQEMGKLETFKRSRIAYSDEVNELLGDAGNSSENQRAHLLDNTERLERSRRLEAGYQIAV
ETEIQIGEMLENLSDREKIQRARDRLRDADANLGKSSRILTGMLRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001293686

ORF Size: 561 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001293686.1](#), [NP_001280615.1](#)

RefSeq Size: 4610 bp

RefSeq ORF: 564 bp

Locus ID: 53611

UniProt ID: [O89116](#)

Cytogenetics: 19 D2

MW: 22.2 kDa

Gene Summary: V-SNARE that mediates vesicle transport pathways through interactions with t-SNAREs on the target membrane. These interactions are proposed to mediate aspects of the specificity of vesicle trafficking and to promote fusion of the lipid bilayers. Involved in vesicular transport from the late endosomes to the trans-Golgi network. Along with VAMP7, involved in a non-conventional RAB1-dependent traffic route to the cell surface used by KCNIP1 and KCND2. May be concerned with increased secretion of cytokines associated with cellular senescence. [UniProtKB/Swiss-Prot Function]