

## Product datasheet for **MR228365**

### Ascl5 (NM\_001270609) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ascl5 (NM\_001270609) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ascl5  
**Synonyms:** Ascl5-ps; bHLHa47; Gm197  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR228365 representing NM\_001270609  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAACAGTAACCTCTGCCGAGCTCTGGTGGACCGTGGGCCCCGGGGGCATGCAGCTAGGTGTCGTGG  
CCCCTGCTGGACAGACGCCGCTGGCAGCCACCGAACCTTGAGCAACGTACCTTTCTATTGTACCCCGG  
CCACTCCGAGCCGCTTACTACGACGCCTACACTGGGGTGTCCCTACGTGCCCTTCCCGGCGCCTTC  
GGGTTTACGACTACCCGTTTCGAGCCTGCCTTCATCCAGAAGCGCAACGAGCGAGAGCGGCAGCGCTCA  
AGTGCGTGAATGAGGGCTACGCGCGCCTCCGCGCCACCTCCCCGCGCCCTGACTGAGAAGCGGCTCAG  
CAAGGTGGAGACTGCGCGCCGCTCCGCTACATCAAGTACCTGCAGGAGCTGTTGAGCGCCACCCCC  
GACGGCGCGCCACCTCCCGCCACCTCCCGCGCGCCGCCACACGGGCCACAGCAACGTGCCGAGCCCT  
CTTCCTTGGTGGCCGAGTCCAGTGGATCCCCCTTTCTCGTCGCTTTTGGAGTCGGAGGAACCCAG  
CCTT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR228365 representing NM\_001270609  
Red=Cloning site Green=Tags(s)

MNSNFCRALVDRGPPGGMQLGVVAPAGQTPLAATEPLSNVPFLLYPGHSEPPYYDAYTGVFPYVPFPGAF  
GVYDYPFEPAFIQKRNERERQVKCVNEGVARLRGHLPGALTEKRLSKVETLRAAIRYIKYLQELLSATP  
DGAPPAATSPPAHTGHSNVPQPSLVAESSGSPFSSSPFLESEEPSL

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM\_001270609

ORF Size: 564 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM\\_001270609.1](#), [NP\\_001257538.1](#)

RefSeq Size: 567 bp

RefSeq ORF: 567 bp

Locus ID: 226439  
Cytogenetics: 1 E4  
MW: 20.8 kDa

## Product images:



Circular map for MR228365