

Product datasheet for MR228400

Cyhr1 (NM_001276322) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cyhr1 (NM_001276322) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cyhr1
Synonyms: 1110031M01Rik; AU042374; Chrp
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228400 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTTTCCAACCAGGACTGAATGGAGCACTGTCTGTCCCACCTGGTATTGGCAGGGGTGTCCCTGC
 ATGCAGCGGTGAGCTCTGTACAGTCCAGTCGAGGAGCAGCTGCCGGTTCCTGCTCCAGACCTTTGCTGC
 CATCATTATGTTGGCCCCAGGGCCGAGCACACATGAAGACTGTCTCGCTGGAGCTGGGTGCCACAGTC
 ATCGGTCTGCCTTTCTGGCCTTCGATTTTCACTGGGTGAATGGGACCGCTTTCTGCCAACCTGCTTC
 TGGGAGGAGGCATGGTGTGGCAGTGGTGGTATCACCTGGGCCCTGAAGGGTGTCTGTCTGCTGGCA
 GGCAGTACTGCTGGTGGTTCAGTGACCATCCTTATTGTGGCTGTTTTACAGCTAACACTTACGGGATG
 TGGGGGGTACCCTGCTGGGTGTGGCAGGCCTTTGAGCAGGCTGGAAGAAGACAGGCTGCTACTGCTGC
 CAAAGGAGGATGTCTGTCTGGTGGGCCCTGGCTGCAGGTAGTTGGGCCATTGCAGGGCCCTGCACACACA
 GCGCCTGCAGTGGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228400 protein sequence
 Red=Cloning site Green=Tags(s)

MVSKPRTEWSTVLSHLVLAGVSLHAAVSSVQSSRGAAAGFLLQTFAAIIMLAPGPSTHEDCLAGAWVATV
 IGLPLLAFDFHWVNGDRSSANLLLGGGMVLAVAGDHLGPEGCSVAGQAVLLVVAVTILIVAVFTANTYGM
 WGGTLLGVAGLLSRLEEDRLLLLPKEDVCRWALAAGSWAYCRALHTQRLQWE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001276322

ORF Size: 579 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001276322.1](#), [NP_001263251.1](#)

RefSeq Size: 1224 bp


RefSeq ORF: 579 bp

Locus ID: 54151

UniProt ID: [Q9QXA1](#)

Cytogenetics: 15 D3
 MW: 20.2 kDa

Product images:



Circular map for MR228400