

Product datasheet for **MR228444**

4833439L19Rik (NM_001252647) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 4833439L19Rik (NM_001252647) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: 4833439L19Rik
Synonyms: 4930558H15Rik; C81457; Kiaa1191; P33monox
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR228444 representing NM_001252647
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGTTCCTGATCACAAAGCAGACCCAGGAGAGCATCCAGCGTTTTGAACAACAGGCAGGACTGAGAG
 ATGCTGGCTACACACCCATAAGGGCCTTACCACGGAGGAGACCAAGTACCTTCGAGTGGCAGAAGCACT
 CCATAAACTAAAGCTGCAGAGTGGAGAGACAGCAAAAGAAGAGAAGCACCCCTGCATCAGCCAGTCCACC
 CCGAGCAGCACCCCTCACGCGTCCCCTAAGCAGAAGTCCAGAGGCTGGTCCCATCTGGCTCCTCCACAG
 CCTTACCTGCTCCGAATCCTCACACCATGGATCCTGGGAGTGGAAATGATAGGAAGTCCAGAGACAAATG
 GAGCCTCTTTGGACCACGACCCCTGCAGAAATCTGATTCTGGTTTGTATTTCAGGCTTACAAAGGTGCC
 CCAAGACCCCTCTCCTATGGAGGTGATGCGTGCCAGGCCACACGAGTAGGTGAAGATCCAGCAACTTTCA
 AACCACCAAAGATGGATGTCCCAATGGTAGAAGGAAAGAAGCAACCACTTCGAACCCATAATCTCAAACC
 CCGTGACCTGAATGTGCTCACACCCACTGGCTTC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228444 representing NM_001252647
 Red=Cloning site Green=Tags(s)


MSSLITKQTQESIQRFEQQAGLRDAGYTPHKGLTTEETKYLRVAEALHKLKLSGETAKEEHPASAQST
 PSSTPHASPKQKSRGWFPSGSSTALPAPNPHTMDPGSGNDRNSADKWSLFGPRPLQKSDSGFAIQAYKGA
 PRPSPMEVMRAQATRVGEDPATFKPPKMDVPMVEGKKQPLRTHNLKPRDLNVLTPTF

SGP**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV


[View online »](#)

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001252647

ORF Size: 594 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001252647.1](#), [NP_001239576.1](#)

RefSeq Size: 2315 bp

RefSeq ORF: 597 bp

Locus ID:	97820
UniProt ID:	Q9DBN4
Cytogenetics:	13 B1
MW:	22.1 kDa
Gene Summary:	Potential NADPH-dependent oxidoreductase. May be involved in the regulation of neuronal survival, differentiation and axonal outgrowth (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228444