

Product datasheet for **MR228476**

Zfp239 (NM_001001792) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Zfp239
Synonyms:	Mok-2; Mok2; Znf239
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR228476 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGCACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCTGAGAAACCTTACAAATGCGACAAGTGCGGGAAGGGTTTACCAGGAGCTCGAGTCTGCTTGTC
ATCATTCGGTCCATACCGGTGAGAAACCTTTCAAGTGTGACAGGTGTGGGAAGGGCTTCAGCCAGAGCTC
AAAGCTCCACATCCACAAGAGAGTCCACACCGGTGAGAAGCCGTATGCGTGCGAGGAGTGCAGTATGAGC
TTCAGTCAGCGATCCAACCTGCACATCCACCAACGTGTCCACACCGGAGAGAGGCCCTACAAGTGTGGGG
AGTGCGGAAAAGGCTTCAGCCAGAGCTCGAACCTCCACATCCACCGGTGCACCCACACGGGAGAGAAGCC
GTACCAAGTGTACGAATGTGGGAAAGGCTTCAGCCAGAGTTCAGACCTTCGGATCCACCTCCGAGTACAC
ACCGGGGAAAAGCCCTACCACTGCGCAAGTGTGGCAGGGCTTCAGCCAGAGCTCCAACTCCTCATCC
ATCAGAGAGTCCACACGGGTGAGAAGCCGTATGAGTGCAGCAAGTGTGGCAAGGGCTTCAGCCAGAGCTC
GAACCTCCACATCCACAGCGGTTACCGCAAGGAGCTTCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228476 protein sequence
Red=Cloning site Green=Tags(s)


MAEKPYKCDKCGKGFTRSSLLVHHSVHTGEKPFKCDRCGKGSQSSKLHIHKRVHTGEKPYACEECGMS
FSQSRNLHIHQRVHTGERPYKCGECGKGSQSSNLHIHRCRTHTEKPYQCYECGKGSQSSDLRIHLRVH
TGEKPYHCGKCGQGFQSSKLLIHRVHTGEKPYECSKCGKGSQSSNLHIHQRVHRKELH

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001001792

ORF Size: 603 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22µm filter is required.

RefSeq:	<u>NM_001001792.2</u>
RefSeq Size:	2828 bp
RefSeq ORF:	606 bp
Locus ID:	22685
UniProt ID:	<u>P24399</u>
Cytogenetics:	6 55.62 cM
MW:	22.8 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228476