

Product datasheet for MR228492

Chtop (NM_001293778) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Chtop (NM_001293778) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Chtop
Synonyms: 2500003M10Rik; Fop; Srag
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228492 representing NM_001293778
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCTGCACAGTCAGCGCCGAAAGTTGTGCTAAAAGCACCACCAAGATGTCTCTAAATGAGCGCTTTA
 CTAATATGCTGAAGAACAACAGCCGATGCCAGTGAATATTCGGGCTTCGATGCAGCAGCAGCAGCAGCT
 AGCCAGTCCAGAAACAGAAGACTGGCCAGCAGATGGAGAATAGACCCTCTGTCCAGGCAGCATTAAAA
 CTTAAGCAGAAGAGCTTAAAGCAGCGCCTGGGTAAAGATAATATCCAGGCACGGTTAGGCCGACCCATAG
 GTGCCCTGGCCAGGGGAGCAATTGGAGGAAGAGGCCTACCATAATCCAGAGAGGCTTGCCCGAGGAGG
 ACTACGTGGGGGACGTGCTACCAGAACCCTGCTTAGGGGTGGGATGTCGCTCCGAGGTCGGGGTATGATA
 GGTCCGGGAAGAGGGGGCTTTGGAGGCAGAGGCCGAGGTCTGGCCGAGGGAGAGGTGCCCTCACTCGCC
 CTGTATTGACCAAGGAGCAGCTGGACAACCAATTGGATGCATACATGTCGAAAACAAAGGACACCTGGA
 TGCTGAATTGGATGCCTACATGGCACAGACAGATCCTGAAACCAATGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR228492 representing NM_001293778
 Red=Cloning site Green=Tags(s)

MAAQSA PKVVLKSTTKMSLNERFTNMLKNKQMPVNI RASMQQQQQLASARNRRLAQOMENRPSVQAALK
 LKQKSLKQRLGKSNIQARLGRPIGALARGAIGRGLPIIQRLPRGGLRGGTRTLRLRGMSLRGRGMI
 GRGRGGFGGRGRGRGRGALTRPVLTKEQLDNQLDAYMSKTKGHLDAELDAYMAQTDPETND

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001293778

ORF Size: 609 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001293778.1](#), [NP_001280707.1](#)

RefSeq Size: 2300 bp

RefSeq ORF: 612 bp

Locus ID: 66511

UniProt ID: [Q9CY57](#)

Cytogenetics: 3 F1

MW: 22.6 kDa

Gene Summary: Plays an important role in the ligand-dependent activation of estrogen receptor target genes (By similarity). May play a role in the silencing of fetal globin genes (PubMed:20688955). Recruits the 5FMC complex to ZNF148, leading to desumoylation of ZNF148 and subsequent transactivation of ZNF148 target genes (PubMed:22872859). Required for the tumorigenicity of glioblastoma cells. Binds to 5-hydroxymethylcytosine (5hmC) and associates with the methylosome complex containing PRMT1, PRMT5, MEP50 and ERH. The CHTOP-methylosome complex associated with 5hmC methylates H4R3 and transactivates genes involved in glioblastomagenesis (PubMed:25284789).[UniProtKB/Swiss-Prot Function]