

## Product datasheet for MR228495

### Plekhb2 (NM\_175421) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Plekhb2 (NM\_175421) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Plekhb2  
**Synonyms:** 2310009M15Rik; AU018239; evt-2; Evt2; Phdc  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR228495 representing NM\_175421  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGTTCGTGAAGAGTGGATGGTTACTTCGGCAGAGCACCATTCTGAAACGCTGGAAGAAGAATTGGT  
TCGACCTGTGGTCAGACGGTCACCTGATCTACTACGATGATCAGACTCGGCAGAGCATAGAGGATAAGGT  
CCACATGCCCGTGGACTGCATCAATATCCGCACGGGCATGAGTGCCGGGACATCCAGCCTCCAGATGGG  
AAGCCCAGAGACTGTCTGCTGCAGATCGTTTGCCGAGACGGGAAGACCATCAGTCTCTGTGCAGAGAGCA  
CAGACGATTGCCTGGCATGGAAGTTTACTGTCAGGATCCAGAACAACACAGCTTACGTTGGTTCAGC  
AATCCTGTCTGAAGAGACTGCAGTGGCCGCGTCCCGCCTCCCTACGCAGCCTATGCTACACCAACCCCT  
GAGGTCTACGGCTATGGTCCATACAGCGGCGCATACCCCGCAGGAACCAAGTTGTCTATGCCGCCAACG  
GGCAGGCATATGCAGTGCATACCAAGTACCCGTATGCAGGTAATGCACGCTGGCTGCCGCTCCCGGTCTC  
TGCGTACAGCGGGGCGAGGAGGAGGAGGCTCCTTTCCAGTTGGAAGTAAA

**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR228495 representing NM\_175421  
Red=Cloning site Green=Tags(s)


MAFVKSGWLLRQSTILKRWKNWFDLWSDGHLIYDDQTRQSIEDKVHMPVDCINIRTGHECRDIQPPDG  
KPRDCLLQIVCRDGTISLCAESTDDCLAWKFTLQDSRTNTAYVGSAILSEETAVAASPPPYAAYATPTP  
EYVYGPYSGAYPAGTQVYYAANGQAYAVPYQYPYAGNARWLPVPSAYSARRRRPPFQLEVK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_175421

**ORF Size:** 612 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_175421.1](#), [NP\\_780630.1](#)

**RefSeq Size:** 2794 bp

**RefSeq ORF:** 615 bp

**Locus ID:** 226971

**Cytogenetics:** 1 B

**MW:** 23.4 kDa

**Gene Summary:** Involved in retrograde transport of recycling endosomes.[UniProtKB/Swiss-Prot Function]