

Product datasheet for MR228519

Clec2f (NM_001277202) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Clec2f (NM_001277202) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Clec2f
Synonyms:	Clec2f-ps; Clr-c; Clrc; Clrz
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR228519 representing NM_001277202 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGACAGCCTCACAGGACGAAGTGGCTCACGAGAGATTGCTGAATACAGAGGTCCTACACAGACCACC
 TAGAAATGGAGAACAGGTAAAAAATCCAAGAAAAATTTGTGGAATTGTCTCTTCCAAGTCTCATTG
 TCGGCGTTTCCGCTGCTGTGGAGTGACCATGGTCTCATTGTAGTTGTGGTCTCATTGTAGTTGTAGTT
 GCTTTGTCAGCAATAATGATAAACAGATGTCAACCAAAATACCTATGCTGCTTGCTGCAAACTGGA
 CTGGAGTTGGAATAAATATTTTATTTTCTGAATACACAAGTAAGTGGACATTTGCCAGGCCTTCTG
 CATGGCCCAAGAGGCCCACTAGCTCGGTTTGACAACGAGGAGGAGCTGAATTTCTAAAGAGACACATG
 AATTCTAGCCACTGGATTGGACTACACAGAGACTCGTCAGAGCACCTTGAGGTGGACAGACAACACTG
 AGTATAACAACACGTTTTTATTCAAGGAGATGGAGAATGTGGCTTCTGAGTGACAATGGGATCAGCAG
 TAGCAGGGACTACATACCTAGGAAGTGGATTGTAGCCGGTCCAGCAACTATATGTTACAGTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>MR228519 representing NM_001277202 Red=Cloning site Green=Tags(s)
-------------------	---

MTASQDEVAHERLLNTEVTTTDHLENGETGKKIQEKFCGIVSSKSHCRRFRCCGVTMVLIVVVVLIVVVV
 ALSAIMIKQMSTKNITYAACLNWTGVGNKYFYFSEYTSNWTFAQAFCAQEAQLARFDNEEELNFLKRHM
 NSSHWIGLHRDSSEHPWRWTDNTEYNNFTLIQGDGECGFLSDNGISSRDYIPRKWICSRSSNYMLQC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001277202

ORF Size: 624 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001277202.1](#), [NP_001264131.1](#)

RefSeq Size: 657 bp

RefSeq ORF: 627 bp

Locus ID: 435921
 Cytogenetics: 6 F3
 MW: 24.3 kDa

Product images:



Circular map for MR228519