

Product datasheet for **MR228532**

Gm11568 (NM_001205030) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gm11568 (NM_001205030) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gm11568
Synonyms:	OTTMUSG00000002207
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR228532 representing NM_001205030 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACCAACTCCTGCTGCTCCCCTTGCTGCCAGCCACCTGCTGCAGGACCACCTGCTGCAGGACCACCT
 GCTGGAGACCAAGCTGTGTGACCAGCTGTTGCCAGCCCTGCTGCCAGCCAGCTGCTGTGGCTCCAGCTG
 CTGCCAGCCTTGCTGTCAAACCACCTGCTGCAGGACCTGCTTCCAGCCATGCTGTGTGAGCAGCTGCTGC
 CGCACACCCTGCTGCCAGCCCTGCTGCTGTGTGCCAGCTGCTGCCAGCCCTGCTGCCAGCCAGCTGTT
 GTCAGTCCAGCTGCTGCCAGCCAGCTGCTGCCAGCCAGCTGCTGTGAGTCCAGCTGCTGCCAGCCAG
 CTGCTGTGAGTCCAGCTGCTGTGAGCCAGGTGCTGTATCTCCAGCTGCTGTGAGCCAGGTGCTGTATC
 TCCAGCTGCTGCCAGCCCTGCTGCAGGCCCTAGCTGCTGTGAGTCCAGCTGCTGCAAGCCCTGCTGCCAGC
 CTTTCTGCCTCAATCTGTGCTGCCAGCCAGCCTGCTCTGGACCTGTGACCTGCACCAGGACTTGCTACCA
 GCCCACCTGCGTCTGTGTGCCTGGCTGCCTGTCCAAGGCTGTGGGTCCAGCTGCTGTGAGCCCTGTGGC
 TGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)


Protein Sequence: >MR228532 representing NM_001205030
 Red=Cloning site Green=Tags(s)

MTNSCCSPCCQPTCCRTTCRTTCWRPSCVTSCCQPCCPSCCGSSCCQPCCTTCRTCFQPCCVSSCC
 RTPCCQPCCVSSCCQPCCPSCCQSSCCQPCCPSCCQSSCCQPCCPSCCQSSCCQPRCCISSCCQPRCCI
 SSCCQPCCRPSCCQSSCCQPCCPFCFLNLCCQACSGPVTCTRTCYQPTCVCPGCLSQGCGSSCCPEPCG
 C

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001205030

ORF Size: 633 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001205030.1, NP_001191959.1

RefSeq Size: 1144 bp


RefSeq ORF: 636 bp

Locus ID: 432600

Cytogenetics: 11

MW: 22.6 kDa

Product images:



Circular map for MR228532