

## Product datasheet for **MR228687**

### Shisa5 (NM\_001284332) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Shisa5 (NM\_001284332) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Shisa5  
**Synonyms:** 2310008D10Rik; 6430628I05Rik; mShisa5; Scotin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR228687 representing NM\_001284332  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCTGCGCCGGCGCCCTCTCTGTGGACCCTATTGCTGTGCTGTTGCTGCTGCCGCGCCTCCGGGTG  
CCCATGGTGAGCTGTGCAGGCCCTTTGGTGAAGACAATTCGATCCAGTGTCTGTCTGATTTCTGTTG  
TGTTCTCTGTTCCAACCAATACTGCTGCTCGGACGTGCTGAGGAAAATCCAGTGGATGAGGAAATGTGT  
CCTGAGCCAGAGTCCAGCAGATTTCCACCCCGCGGAGGAGACACCCGAACATCTGGGTTCCAGCGCTGA  
AATTTTCGATCCAGTTTTGACAGTGACCCTATGTACGGTTTCGGAGCGACCGTCCGATTGGCGTGACCAT  
CTTTGTGGTGTATTGCCACTATCATCTGCTTACCTGCTCCTGCTGCTGTCTGTATAAGATGTGC  
TGCCCCAACGCCCTGTCGTGACCAACACCACAATACTACCGTGGTTCATGCCCTTACCCTCAGCCTC  
AACCTCAACCTGTGGCCCCAGCTATCCTGGACCAACATACCAGGGTACCATCCCATGCCCCCCAGCC  
AGGAATGCCAGCAGCACCTACCCAACGCAGTACCACCACCTACCTGGCCAGCCACAGGGCCGCCA  
CCCTACCATGAGTCCTTGGCTGGAGCCAGCCAGCCTCCATACAACCCGACCTACATGGATTCCCTAAGA  
CAATTCCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR228687 representing NM\_001284332  
 Red=Cloning site Green=Tags(s)

MAAPAPSLWLLLLLLLLPPPPGAHGLCRPFGEDNSIPVFCPDFCCGSCSNQYCCSDVLRKIQWNEEMC  
 PEPSSRFSTPAEETPEHLGSAKFRSSFDSDPMSGFGATVAIGVTIFVVFIIATIIICFTCSCLYKMC  
 CPQRPVVTNTTTTTVVHAPYPQPQPVPVAPSYGPTYQGYHPMPQPQMPAAPYPTQYPPPYLAQPTGPP  
 PYHESLAGASQPPYNPTYMDSLKTIP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001284332  
 ORF Size: 708 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001284332.1</a> , <a href="#">NP_001271261.1</a>
<b>RefSeq Size:</b>	1923 bp
<b>RefSeq ORF:</b>	711 bp
<b>Locus ID:</b>	66940
<b>UniProt ID:</b>	<a href="#">Q9D7I0</a>
<b>Cytogenetics:</b>	9 F2
<b>MW:</b>	26.1 kDa
<b>Gene Summary:</b>	Can induce apoptosis in a caspase-dependent manner and plays a role in p53/TP53-dependent apoptosis.[UniProtKB/Swiss-Prot Function]