

## Product datasheet for MR228737

### Rgs19 (NM\_001291205) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Rgs19 (NM_001291205) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rgs19
Synonyms:	2610042F04Rik; AI324841; AW547781; GAIP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR228737 representing NM_001291205 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTATCTGCCAGGATCCTATGCTGCAGTGAGGGTCCAGAGCTGGCACCACAGAATCTGGCAACCCCC  
TCCACTACCCATGCCCACCCACATGAGGCTGAGAAACAGCACACAGGGCCAGAGGAGGCAGACAGGCC  
TCCCTCGATGTCCAGCCATGACGCAGCCCCCTCAGGACCTCCAGTCGCAATCCCTGCTGCTTATGCTGG  
TGCTGCTGCTGTAGCTGCTCCTGGAACCAAGAACGGCAACGAGCTTGGCAGGTTTCTCGGAAAGCAAGC  
TGCAACCCCTCCCAGCTGTGAAGTCTGCACTCCACCAAGCCCTAAGGAAGTACAGAGCTGGGCACAGTC  
CTTTGACAAGTTGATGCATAGCCCCACGGGCCGAGTGTATTCCGGGCATTCTGCGCACAGAATACAGC  
GAAGAGAACATGCTCTTCTGGCTGGCTTGTGAGGAGTTGAAAGCAGAGGCCAACCAACATGTGGTGGACG  
AGAAGGCGCGACTTATCTATGAGGACTACGTGTCCATCCTGTCCCCAAGGAGGTAAGCCTGGACTCCCG  
TGTGCGAGAAGGCATCAATAGGAAGATGCAAGAGCCATCGCCACACACATTTGATGATGCACAGCTGCAG  
ATCTACACCCTCATGCACCGGACTCGTACCCTCGATTCTTACCTCCCCACCTACCGCTCTCTATTAC  
TCCAGGGGGCCCCACAGTCTCTGAGGCC

**ACGCGT**ACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR228737 representing NM\_001291205  
 Red=Cloning site Green=Tags(s)

MVSARILCCSEGPPELAPQNLGNPLHSPMPTPHEAEKQHTGPEEADRPPSMSSHDAAPSGPPSRNPCCLCW  
 CCCCSCSWNQERQRAWQVSRESKQLPLPSCVCTPPSPKEVQSWAQSFDKLMHSPTGRSVFRAFLRTEYS  
 EENMLFWLACEELKAEANQHVVDEKARLIYEDYVSILSPKEVSLDSRVREGINRKMQEPSPHTFDDAQLQ  
 IYTLMHRDSYPRFLTSPTYRSLLLQGAPQSSEA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001291205  
 ORF Size: 729 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001291205.1</a> , <a href="#">NP_001278134.1</a>
<b>RefSeq Size:</b>	1512 bp
<b>RefSeq ORF:</b>	732 bp
<b>Locus ID:</b>	56470
<b>Cytogenetics:</b>	2 103.72 cM
<b>MW:</b>	27.9 kDa
<b>Gene Summary:</b>	Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form. Binds to G-alpha subfamily 1 members, with the order G(i)a3 > G(i)a1 > G(o)a >> G(z)a/G(i)a2. Activity on G(z)-alpha is inhibited by phosphorylation and palmitoylation of the G-protein (By similarity). [UniProtKB/Swiss-Prot Function]