

## Product datasheet for **MR228779**

### **Bdnf (NM\_001285418) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Bdnf (NM\_001285418) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Bdnf  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR228779 representing NM\_001285418  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACCATCCTTTTCCTTACTATGGTTATTTCACTTCGGTTGCATGAAGGCGGCCCCATGAAGAAG  
TAAACGTCCACGGACAAGGCAACTTGGCCTACCCAGGTGTGCGGACCCATGGACTCTGGAGAGCGTGAA  
TGGGCCAGGGCAGGTTTCGAGAGGTCTGACGACGACATCACTGGCTGACACTTTTGAGCAGTCATCGAA  
GAGCTGCTGGATGAGGACCAGAAGGTTTCGCCCCAACGAAGAAAACCATAAGGACGCGGACTTGTACTT  
CCCGGGTGATGCTCAGCAGTCAAGTGCCTTTGGAGCCTCCTACTCTTTCTGCTGGAGGAATACAAAA  
TTACCTGGATGCCGAAACATGTCTATGAGGGTTCGGCGCCACTCCGACCCTGCCGCCGTGGGAGCTG  
AGCGTGTGTGACAGTATTAGCGAGTGGGTACAGCGGCAGATAAAAAGACTGCAGTGGACATGTCTGGC  
GGACGGTCACAGTCTAGAGAAAGTCCCGGTATCCAAAGGCCAACTGAAGCAGTATTTCTACGAGACCAA  
GTGTAATCCCATGGGTTACACCAAGGAAGGCTGCAGGGGCATAGACAAAAGGCACTGGAATCGCAATGC  
CGAACTACCAATCGTATGTTCCGGCCCTTACTATGGATAGCAAAAAGAGAATTGGCTGGCGATTATAA  
GGATAGACTTCTGTGTATGTACTGACCATTAAGGGGAAGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR228779 representing NM\_001285418  
**Red=Cloning site Green=Tags(s)**

MTILFLTMVISYFGCMKAAPMKEVNVHGQGNLAYPGVRTHGTLESVNGPRAGSRGLTTTSLADTFEHVIE  
ELLDEDQKVRPNEENHKDADLYTSRVMSSQVPLEPPLFLLEEYKNYLDAAANMSMRVRRHSDPARRGEL  
SVCDISSEWVTAADKKTAVDMSGGTVTVLEKVPVSKGQLKQFYETKCNPMGYTKEGCRGIDKRHWNSQC  
RTTQSYVRALTMSKRIGWRFIRIDTSCVCTLTIKGRG

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)



<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001285418.1</a></u> , <u><a href="#">NP_001272347.1</a></u>
<b>RefSeq Size:</b>	3866 bp
<b>RefSeq ORF:</b>	750 bp
<b>Locus ID:</b>	12064
<b>UniProt ID:</b>	<u><a href="#">P21237</a></u>
<b>Cytogenetics:</b>	2 56.63 cM
<b>MW:</b>	28.6 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the nerve growth factor family. It is involved in the growth, differentiation and survival of specific types of developing neurons both in the central nervous system (CNS) and the peripheral nervous system. It is also involved in regulating synaptic plasticity in the CNS. Expression of a similar gene in human is reduced in both Alzheimer's and Huntington disease patients. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar processing to generate mature protein. [provided by RefSeq, Oct 2015]