

## Product datasheet for **MR228806**

### Adal (NM\_001290812) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Adal (NM\_001290812) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Adal  
**Synonyms:** 4930578F03Rik  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR228806 representing NM\_001290812  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAATTGACAGAAGAGGTGGCCCAACAATAGCCAGGGAGACTGTGCGAACTTGCTAAAGAGTTCTTCC  
TCTCTACTGAGAATACAGTTCTTGGCCTTGACCTCAGTGGAGATCCAACATAGGACAAGCTAACGACTT  
CTTGGAGCCTCTTTTAGAAGCTAAAAAGCAGGTCTGAAGTTGGCATTGCATCTTGCCGAGATCCAAAC  
AGGGAAAAAGAAAATCAGATGCTGCTGAGTCTGCTCCCTGACAGAATTGGGCATGGGACATTCCTTAGCG  
CCTCTGAGGCAGGGCCCTGGACCAAGTGGACTTCGTGAGGCAACACCAGATACCCCTCGTGTGGAAAGC  
AAATGCAGAGCTCTACATGTGCAGGGCAGGCCCTCTTCTGCTGAGGCACATCCCCACTCAGAACAAAGTCC  
ACATCACAAGTAATGAACATTTTGTGGAACTTTGTTTGACCTCAAACATCAAAGCCAGACAGTCCCTT  
CTTATGATCAGCATCACTTTGGATTTTGGTACAGCATTGCTCATCCTTCTGTGATCTGTACTGATGATAA  
AGGTGTTTTTGGACATACCTTTCTCAAGAATATCAACTGGCAGCAGAAACATTTAATTTGACCCCATTT  
CAGGTGTGGGATCTGTCCTATGAATCCATCAACTACATCTTTGCTTGTGACAACACCAGATCTGAACTGA  
GAAAAAGGTGGACTCACTGAAGCAGAAAGTGTAAACTGTAATGAGGTGAACTACTTT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR228806 representing NM\_001290812  
Red=Cloning site Green=Tags(s)

MAIDRRGGPTIARETVELAKEFFLSTENTVLGLDLSGDPTIGQANDFLEPLLEAKKAGLKLALHLAEIPN  
 REKENQMLLSLLPDRIGHGTFLLSASEAGALDQVDFVRQHQIPLVLEANAELMCRAGPLLLRHIPQTQKS  
 TSQVMNILWELCLTSNIKSQTVPSYDQHFFGFWYSIAHPSVICTDDKGVFATYLSQEYQLAAETFNLTPF  
 QVWDLSEYSEINYIFACDNTRSELKRWTHLKQKVLNCNEVNYF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001290812  
**ORF Size:** 759 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001290812.1</a> , <a href="#">NP_001277741.1</a>
<b>RefSeq Size:</b>	2369 bp
<b>RefSeq ORF:</b>	762 bp
<b>Locus ID:</b>	75894
<b>UniProt ID:</b>	<a href="#">Q80SY6</a>
<b>Cytogenetics:</b>	2 E5
<b>MW:</b>	29.1 kDa
<b>Gene Summary:</b>	Putative nucleoside deaminase. May catalyze the hydrolytic deamination of adenosine or some similar substrate and play a role in purine metabolism (By similarity).[UniProtKB/Swiss-Prot Function]