

Product datasheet for MR228822

Folr1 (NM_001252554) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Folr1 (NM_001252554) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Folr1
Synonyms:	FBP1; Folbp-1; Folbp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR228822 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGGCTCACCTGATGACTGTGCAGTTGTTGCTCCTGGTATGTGGATGGCCGAATGTGCTCAGTCCAGAGCTACTCGGGCCAGGACTGAACTTCTCAATGTCTGCATGGATGCCAAACACCACAAAGAAAACCGGGCCC TGAGGACAATTTACACGACCAGTGCAGCCCTGGAAGACGAATTCCTGCTGTTCCACGAACACAAGCCAG GAAGCACATAAGGACATTTCTACCTGTACCGTTCAACTGGAACCACTGCGGAATATGACATCGGAAT GCAAACGGCACTTTATCCAAGACACCTGCCTCTATGAGTGTCCCGAATTTGGGACCCTGGATCCAGCA GGTGGACCAGAGCTGGCGCAAAGAGCGGATCCTTGATGTCCCTGTGCAAAGAGGACTGTGAGCAGTGG TGGGAGGACTGCCAGAGCTCTTTTACCTGCAAGAGCAATTGGCACAAGGGATGGAATGGTCTCGGGGC ATAACGAGTGTCTGTGGGAGCCTCCTGCCATCCCTTACCTTCTACTTCCCACATCTGCTGCTGTGTG TGAGGAAATCTGGAGTCACTCCTACAAGCTCAGCAACTACAGTCGAGGGAGCGGCCGCTGCATTCAGATG TGGTTCGACCCAGCCAGGGCAACCCCAACGAGGAAGTGGCGAGGTTCTATGCCGAGGCCATGAGTGGAG CTGGGTTTCATGGGACCTGGCCACTCTGTGCAGCCTGTCTTAGTGCTGCTCTGGGTGATCAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR228822 protein sequence
 Red=Cloning site Green=Tags(s)



MAHLMTVQLLLL VMWMAECAQSRATRARTELLNVCMDAKHHKEKPGPEDNLHDQCSPWKTNSCCSTNTSQ
 EAHKDISYLYRFNWNHCGTMTSECKRHF IQDTCLYECSPNLGPWIIQQVDQSWRKERILDVPLCKEDCQQW
 WEDCQSSF TCKSNWHKGNWSSGHNECPVGASCHPFTFYFPTSAALCEEIWSHSYKLSNYSRGSGRCIQM
 WFDPAQGNPNEEVARFYAEAMSGAGFHGTWPLLCSLSLVLLWVIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001252554

ORF Size: 765 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001252554.1](#), [NP_001239483.1](#)

RefSeq Size: 1324 bp

RefSeq ORF: 768 bp

Locus ID: 14275


UniProt ID: [P35846](#)

Cytogenetics: 7 E2

MW: 29.4 kDa

Gene Summary: Binds to folate and reduced folic acid derivatives and mediates delivery of 5-methyltetrahydrofolate and folate analogs into the interior of cells. Has high affinity for folate and folic acid analogs at neutral pH. Exposure to slightly acidic pH after receptor endocytosis triggers a conformation change that strongly reduces its affinity for folates and mediates their release. Required for normal embryonic development and normal cell proliferation. Required for renal folate reabsorption.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR228822